








Brand Name of Product	LumCheck™
Generic Name of Product	Lumen cleaning challenge device.
Product Code Number(s)	WLC-101, WLC-102
Purpose of Product	Challenge automated lumen cleaning function of cleaning machines.
Range of Applications for Product	May be used in cleaning equipment intended to clean lumened instruments.
Key Specifications of Product	<ul style="list-style-type: none"> • Coupon with a test soil comprised of simulated blood soil. • Single-channel lumen holder for a lumened instrument. The holder unscrews for placement of the test coupon.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	Store in the protective foil pouch at a temperature of 2- to 25 °C (35.6- to 77 °F).
Packaging Contents	<ul style="list-style-type: none"> • Each coupon is placed in an individual water-tight foil pouch. • Twenty-five (25) tests per box. • Holder sold separately.
Shelf Life	Eighteen (18) months from date of manufacture (Mfr.). (Note: Consult package for expiration date.)

Instructions for Using Product	
Description of Use(s)	Challenge automated lumen cleaning function of cleaning machines.
Preparation for Use	N/A
Diagrams (drawings, pictures):	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Unscrew the LumCheck™ holder. 2. Detach part A (marked with “slot”) from part B. 3. Open the protective pouch of LumCheck™ and insert the test coupon (part C) in part B. Do not touch the area covered with the test soil. 4. Close LumCheck™ holder by screwing part A into part B. 5. Connect the holder with one of the channel irrigation systems of the washer (luer-lock). <i>(Note: Consider using a minimum of one (1)-LumChek™ test per level if utilizing a multi-level ultrasonic cleaner).</i> 6. Start the wash cycle according to Mfr.’s IFU. 7. At the end of the cycle, retrieve the holder from the washer. 8. Open the holder by unscrewing part A from part B and remove the test coupon without touching the soiled areas. 9. Visually inspect the coupon for residual soil. Compare the test result to the interpretation chart, below.

Interpretation of Results							
	Result	Optimum Results No Residue Remaining	Visible Fibrin	Moderate/High Level of Fibrin	Visible Hemoglobin	Visible Hemoglobin and Visible Fibrin	Moderate/High Level of Fibrin and Hemoglobin
	Parameters that may affect LumCheck Results:		Detergent Selection and/or Dosing	Detergent Selection and/or Dosing	Water Temperature and/or Flow	Detergent Selection and/or Dosing Water Temperature and/or Flow	Detergent Selection and/or Dosing Water Temperature and/or Flow
Contraindications of Test Results	N/A						
Documentation	Record results on a log sheet .						
Special Warnings and Cautions	N/A						
Disposal	Dispose of in a ☒ biohazard container.						

Reprocessing Instructions	
Point of Use	N/A
Preparation for decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	TOSI®
Other Product Support Documents	N/A
Reference Documents	N/A
Customer Service Contact	Healthmark, A Getinge company 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com