






**Instructions for Use: All-in-One Microscope Camera**

<b>Brand Name of Product</b>	iLoupe XL2 All-in-One Microscope Camera
<b>Generic Name of Product</b>	Microscope Camera
<b>Product Code Number(s)</b>	26700-520
<b>Purpose of Product</b>	Designed for professional optical inspection of items and other equipment.
<b>Range of Applications for Product</b>	
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• One of the lightest, most powerful, and compact mobile microscopes in the world</li> <li>• Palm-sized microscope that does not require any additional monitors to view subjects, giving freedom to take it anywhere you need magnification</li> <li>• Full color display allows for the viewing of images at high magnification, the ability to take photos, the ability to record videos, and the ability to display output on a PC</li> <li>• May also be used as a web camera</li> <li>• Four bright, built-in LEDs with intensity control allows for control of the illumination</li> <li>• The included software has powerful measurement function to measure length, perimeter, angle, and area of the object in view</li> <li>• 2" LCD video screen</li> <li>• Wide magnification range of 3X – 300X</li> <li>• Built-in 5 mega-pixel camera</li> <li>• Features: 180° rotated lens tube, SD card for storing and transferring images; easy-to-use touchpad controls, rechargeable lithium-ion battery, TV video output</li> </ul>

**Shipping & Storage**

<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	
<b>Shelf Life</b>	

**Instructions for Using Product**

<b>Description of Use (s)</b>	Use for inspection of items and other equipment.
<b>Preparation for Use</b>	
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Press the power button located in the upper right-hand corner of the microscope until the screen lights up and images appear.</li> <li>2. Use the high and low magnification control on the left-hand side of the microscope to adjust to achieve desired magnification/focus.</li> <li>3. To take a snapshot of the image on the screen, press the snapshot button located on the left-hand side of the top of the microscope. An hourglass figure will appear and then your snapshot will appear on the screen.</li> <li>4. To view photos, press the  button on the lower right-hand side of the microscope. Navigate through the photos using the left and right arrows of the touchpad.</li> <li>5. To exit the saved photos, pressure the  button again.</li> <li>6. The function menu can be accessed by pressing the home/trash can button. Navigate through the function menu with the arrows on the touchpad. Select the desired function by pressing the center of the touchpad.</li> </ol>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	

<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	
<b>Storage</b>	
<b>Additional Information</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	ProSys™ Instrument Care Product Brochure, ProSys™ Instrument Care Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>