


<b>Brand Name of Product</b>	Instrument Protector Card												
<b>Generic Name of Product</b>	Instrument Protector Card												
<b>Product Code Number(s)</b>	IPC001, IPC002, IPC003												
<b>Intended Use</b>	To protect delicate items.												
<b>Range of Applications for Product</b>	Various sized items.												
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>Made from solid bleached sulfate (SBS)– 265 g/m<sup>2</sup>.</li> <li><i>Paper:</i> Expandable base accommodates curved/larger instruments.</li> <li>Single use.</li> </ul> <p>Instrument Protector Sizes:</p> <table border="1"> <thead> <tr> <th>Instrument Name</th> <th>Instrument Size</th> <th>Dimensions</th> </tr> </thead> <tbody> <tr> <td>IPC001</td> <td>Large</td> <td>5.5- x 9.50 inches</td> </tr> <tr> <td>IPC002</td> <td>Medium</td> <td>3.5- x 6.75 inches</td> </tr> <tr> <td>IPC003</td> <td>Small</td> <td>1.75- x 8.25 inches</td> </tr> </tbody> </table>	Instrument Name	Instrument Size	Dimensions	IPC001	Large	5.5- x 9.50 inches	IPC002	Medium	3.5- x 6.75 inches	IPC003	Small	1.75- x 8.25 inches
Instrument Name	Instrument Size	Dimensions											
IPC001	Large	5.5- x 9.50 inches											
IPC002	Medium	3.5- x 6.75 inches											
IPC003	Small	1.75- x 8.25 inches											

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	<ul style="list-style-type: none"> <li>Store in dry conditions (20 – 60 rh%) at indoor room temperature 15- to 30 °C (59- to 86 °F).</li> <li>Store in original packaging or in an equivalent clean barrier to maintain cleanliness.</li> </ul>
<b>Packaging Contents</b>	100 cards each box.
<b>Shelf Life</b>	N/A

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	For use in sterile barrier systems; to protect items during sterilization, storage, and transportation.
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures)</b>	 <p><b>Figure 1</b></p> <p><b>Figure 2</b></p>




**Figure 3a**





**Figure 3b**



**Figure 3c**



**Figure 3d**



**Figure 3e**



**Figure 4**

**Steps for Use of Product**

1. Fold bottom flap 180°. (**Fig 1**).
2. Fold the locking tabs and insert into the two-tab receiving slots at the bottom of the Instrument Protector. (**Fig 2**).
3. Place instrument into the Instrument Protector.
  - a. Feed instrument through large “top hat” slot to the backside of card. (**Fig. 3a**).
  - b. Next, feed instrument upward through the small “top hat” slot to the front of the card. (**Fig. 3b**).
  - c. Position instrument by working the distal end of the instrument into bottom fold-over pocket. (**Fig. 3c**).
  - d. Open instrument (if applicable) so that its arms extend to outer most boundary of the large “top hat” slot. (**Fig. 3d**).

	<p>e. Fold tongue down over instrument arms and through the “top hat” slot (pull tight to secure instrument). <b>(Fig. 3e).</b></p> <p>4. If using a peel pack, place Instrument Protector Card with finger hole oriented toward the chevron side. <b>(Fig. 4).</b></p>
<b>Interpretation of Test Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	Instrument Protector Cards are not a barrier to sterilant or a sterile barrier.
<b>Disposal</b>	<i>Single use:</i> Dispose of in compliance with facility guidelines.

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>• <i>Compatible with steam sterilization:</i> <ul style="list-style-type: none"> <li>○ 132 °C (270 °F) for 4 minutes</li> <li>○ 135 °C (275 °F) for 3 minutes.</li> </ul> </li> <li>• <i>Extended cycle:</i> 135 °C for 30 minutes.</li> <li>• Material compatibility with EO sterilization.</li> <li>• Validation testing was performed for single and double pouching.</li> </ul>
<b>Storage</b>	<ul style="list-style-type: none"> <li>• Keep in dry conditions (20–60 rh%) and indoor room temperatures 15- to 30 °C (59- to 86 °F).</li> <li>• Keep in original packaging or in an equivalent clean barrier to maintain cleanliness.</li> </ul>
<b>Additional Information</b>	N/A
<b>Related Healthmark Products</b>	Window Style Tip Protectors, Vented Instrument Guards, Duo-Vented Instrument Guards, Cylindrical Style, Instrument Guards.
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark, A Getinge Company 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com