



Brand Name of Product	ScopeDry Check™
Generic Name of Product	ScopeDryCheck™
Product Code Number(s)	SDC-001
Intended Use	To test lumens/channels for detectable levels of residual moisture.
Range of Applications for Product	For any lumened/channeled items.
Key Specifications of Product	<ul style="list-style-type: none"> • Measuring range: Detection down to 50 nL (or 0.05µL) of residual moisture. • ScopeDryCheck™ Card - 3.5- (W) x 2 inches (H). • Response time: 2 to 5 seconds.

Shipping & Storage	
Shipping Conditions & Requirements	Avoid exposure to moisture.
Storage Conditions	<ul style="list-style-type: none"> • Store in a cool dry place at 60-80°F with 30-60% rh. • Avoid exposure to moisture. • Avoid exposure to UV light.
Packaging Contents	N/A
Shelf Life	Two (2) years from the date of manufacture.

Instructions for Using Product	
Description of Use(s)	To test lumens/channels for detectable levels of residual moisture.
Preparation for Use	N/A
Diagrams (drawings, pictures)	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>Figure 1</p> </div> <div style="text-align: center;">  <p>Figure 2A</p> </div> </div>




Figure 2B




Figure 2C







Figure 3

	  
<p>Steps for Use of Product</p>	<ol style="list-style-type: none"> 1. With dry gloved hands, after drying any lumen/channeled items according to the manufacturer’s IFU, with the same air source used to dry the scope, flush air through the biopsy port. 2. On the flipside of the card, with a gloved hand, hold the ScopeDryCheck™ in one of the corners with the blue side up, ready to test if the lumen/channel is dry, while blowing air, aim the scope one (1)-cm away from the card in the following sequence: <ol style="list-style-type: none"> A. <i>Distal end</i>: Flush air through the biopsy port. (Fig. 1, 2A). B. <i>Suction connector</i>: Flush air through the suction valve port. (Fig. 2B). C. <i>Auxiliary water inlet</i>: Flush air through the suction valve port. (Fig. 2C). 3. If any white dots are visible at any time during this test procedure: <ol style="list-style-type: none"> a. resume drying of the item. b. Retest with a new card. (Fig. 3). 4. Repeat test until no white dots are visible.
<p>Interpretation of Test Results</p>	<p>If white dots are visible on the card, the item is not dry. Repeat drying steps and retest with another area of the card to check for residual moisture.</p>
<p>Contraindications of Test Results</p>	<p>N/A</p>
<p>Documentation</p>	<p>N/A</p>
<p>Special Warnings and Cautions</p>	<p>Prolonged exposure to UV light causes the indicator chemical to activate and fade over time.</p>
<p>Disposal</p>	<ul style="list-style-type: none"> • Single use only. • Dispose of according to healthcare facility policy.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	<ul style="list-style-type: none"> • Store in a cool dry place at 60 to 80 °F with 30–60% rh. • Avoid exposure to moisture. • Avoid exposure to UV light.
Additional Information	N/A
Related Healthmark Products	N/A
Other Product Support Documents	ProFormance™ Brochure, ProFormance™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark, A Getinge company 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com