




<b>Brand Name of Product</b>	Gorilla Bags™
<b>Generic Name of Product</b>	Dust Covers
<b>Product Code Number(s)</b>	GB-001 PP, GB-002 PP, GB-003 PP
<b>Intended Use</b>	To be used as protective packaging (as defined by ISO 11607-1)
<b>Range of Applications for Product</b>	Wrapped trays
<b>Key Specifications of Product</b>	Manufactured from spun polypropylene.

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	Clean, dust-free area, and away from fluorescent or ultraviolet light.
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	Five (5)-year shelf life prior to use when stored at: <ul style="list-style-type: none"> <li>• Humidity: 20–80% rh</li> <li>• Temperature: 10 °C (50 °F) to 35 °C (95 °F).</li> </ul>

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	Added protection from environmental conditions as protective packaging.
<b>Preparation for Use</b>	Select proper size bag for wrapped tray.
<b>Diagrams (drawings, pictures)</b>	<div style="display: flex; justify-content: space-around;">   </div> <div style="display: flex; justify-content: center; margin-top: 20px;">  </div>
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Wrap tray (according to policy) with facility approved wrap.</li> <li>2. Once the tray is wrapped, place bag on table.</li> <li>3. Place wrapped tray into Gorilla Bag™. (<b>Fig. 1</b>).</li> <li>4. Fold the opening of the bag forward three times. (<b>Fig. 2</b>).</li> <li>5. Seal with two strips of autoclave tape <ol style="list-style-type: none"> <li>a. Place each strip on the bag</li> <li>b. Wrap tape underneath the fold to secure the closure of the bag. (<b>Fig. 3</b>).</li> </ol> </li> <li>6. After cool down, inspect the Gorilla Bag™ before it is removed from the rack.</li> <li>7. Follow facility polic(ies) with proper handling techniques by lifting the wrapped tray off the rack and set onto the transportation cart to be stored. (<b>Note: Do not drag the tray.</b>)</li> <li>8. After the set has been transported from storage to the OR, remove Gorilla Bag™.</li> </ol>

<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• Is not a microbiological layer.</li> <li>• Never grab/hold the end of the Gorilla Bags™ while transporting/carrying items. This action can compromise the integrity of the bag.</li> <li>• Proper care and handling of the Gorilla bags™ should be observed at all times.</li> </ul>
<b>Disposal</b>	<ul style="list-style-type: none"> <li>• Single-use only.</li> <li>• Landfill or incinerate based upon state or local regulations.</li> <li>• Gorilla Bags™ are composed of polypropylene plastic, which has a recycling code 5 ♻.</li> </ul>

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	Gorilla Bags™ are completely porous to the sterilant and can be used in pre-vacuum steam sterilizers for a maximum temperature of 138 °C (280.4 °F) for a maximum time of 30 minutes.
<b>Storage</b>	N/A
<b>Additional Information</b>	Available color: purple

<b>Related Healthmark Products</b>	Sterilization Packaging Brochure, Sterilization Packaging Price List
<b>Other Product Support Documents</b>	N/A
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark, A Getinge Company 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com