

# XEN

Chemistries  
for the ultimate  
peace of mind

## Rinse



Cart Washer



Washer/  
Disinfector

## Instructions for Use

**XEN Rinse** is a liquid rinse aid for use in automated washer applications, including washer-disinfectors and cart washers. Rinse is low foaming and does not interfere with steam or ethylene oxide sterilization.

Rinse aids the sheeting of rinse water over the cleaned instruments, ensuring the water drains evenly from the surfaces, leaving a spot-free finish and reducing drying time.

Order no.	Description
6036002101	Rinse 1 x 5 L.
6036002201	Rinse 1 x 5 gal.
6036002301	Rinse 1 x 15 gal.



### Directions for Use in Automated Washing Equipment

1. Rinse should be automatically dosed into the final rinse cycle. Dosing is 0.03–0.13 oz./gal. (0.25–1.0 ml/L) at temperatures above 122°F (50°C).
2. Do not use on polycarbonate. Check polysulphone (PSU, PPSU) and similar materials for compatibility before use.

GETINGE 

## Caution

- Do not mix with other chemicals.
- Wear correct PPE and follow applicable guidelines. Additional information on SDS.
- Check material compatibility before use.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

**Danger:** Causes serious eye damage. Causes skin irritation. Wear protective gloves/protective clothing/eye protection/face protection. Wash hands and skin thoroughly after handling. Immediately call a poison center or doctor. For specific treatment, refer to SDS Section 4: First aid measures.

**First Aid:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF ON SKIN: Wash with plenty of water. IF SKIN IRRITATION OCCURS: Get medical advice. Remove contaminated clothing and wash before reuse.

For ingredient information, visit [www.getinge.com](http://www.getinge.com)

## Storage Information

Store in a cool, dry place protected from frost and away from alkalis and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Rinse	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instruments/device manufacturer

Teflon™ is a trademark of Chemours.

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2024 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

**Manufacturer:** Quadralene Ltd, Unit 3, St. Modwen Park, Andresseway Way, Derby, DE21 6YH, United Kingdom.  
**Sales Office, US:** 1 Geoffrey Way, Wayne, NJ 07470, USA  
**Sales Office, Canada:** 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

[www.getinge.com](http://www.getinge.com)

