

XEN

Chemistries
for the ultimate
peace of mind

Chamber Shine+



Cart Washer



Maintenance



Washer/
Disinfectant

Instructions for Use

XEN Chamber Shine+ is a highly concentrated, low-foaming liquid detergent for cleaning and descaling stainless-steel medical instruments and wash chambers using automated washers, including cart washers.



| Order no. | Description |
|------------|----------------------------|
| 6036004201 | Chamber Shine+ 1 x 5 L. |
| 6036004301 | Chamber Shine+ 1 x 5 gal. |
| 6036004401 | Chamber Shine+ 1 x 15 gal. |

Directions for Use in Automated Washing Equipment

1. Chamber Shine+ should be automatically dosed into the wash cycle . Dosing is 0.06-0.26 oz./gal. (0.5–2.0 ml/L), depending on soil load and water quality. Higher doses will be required for heavier scale deposits and for poor water quality.
2. . If staining persists, repeat the process . Ensure thorough rinsing.

GETINGE 

Caution

- Do not mix with other chemicals.
- Wear correct PPE and follow applicable guidelines. Additional information on SDS.
- Check material compatibility before use.

Danger: Causes severe skin burns and eye damage. May be corrosive to some metals. Keep only in original container. Do not breathe fumes, gas, mist, vapors, or spray. Wash hands and skin thoroughly after handling. Wear protective gloves, protective clothing, eye protection, and face protection.

First Aid: IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF INHALED: If breathing is difficult, remove person to fresh air and keep at rest in a comfortable position for breathing. Immediately call a poison center or doctor. For specific treatment refer to section 4 of SDS. Wash contaminated clothing before reuse. IF ON SKIN (or

hair): Immediately remove all contaminated clothing. Rinse SKIN with water or shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison control center or doctor. For specific treatment, refer to section 4 of SDS. Absorb spillage to prevent material damage. Store locked up. Store in corrosive resistant container with a resistant inner liner.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

Contains phosphoric acid and nitric acid. For more ingredient information, visit www.getinge.com

Storage Information

Store in a cool, dry place protected from frost and away from alkalis and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C.

| | Aluminum | Anodized Aluminum | Titanium | Copper | Brass | Gold | Nickel | Cobalt Chrome | Matte chrome plated | Stainless steel | Glass | Ceramics | Rubber | Plastic | Teflon™/ PTFE | Polycarbonate |
|----------------|----------|-------------------|----------|--------|-------|------|--------|---------------|---------------------|-----------------|-------|----------|--------|---------|---------------|---------------|
| Chamber Shine+ | | | | | | | | | | • | • | | | | | |

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instruments/device manufacturer

Teflon™ is a trademark of Chemours.

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2024 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andresseway Way, Derby, DE21 6YH, United Kingdom.
Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA
Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com



MCV00110914 REVA