

## Steriking® LT-Blueline Pouches with Tyvek

*Steriking® Heat sealable pouches and rolls are intended for use as packaging material for reusable instruments in STERRAD® low temperature sterilization by hydrogen peroxide in health care establishments. The products are for single use only. Products are designed to be used by trained healthcare professionals. Benefit of pouches is to allow sterilization, maintain sterility and enable aseptic presentation of packed medical device.*

Preformed sterile barrier systems made of porous material and laminate construction. The pouches are constructed from Tyvek®/plastic films. The heat-sealed pouches are heat sealed prior to processing in the STERRAD® Sterilization Systems.

The device is placed inside the pouch and the four side is sealed.

### Packaging

When measuring the size, it is important to allow extra space for sealing the package and for a flap with which the package can then be easily opened. It is recommended that the packages should be filled to no more than 3/4 of their length. Empty space should be allowed around each instrument. The



medical device should be oriented to ensure aseptic presentation - in other words positioned correctly inside the pack to enable easy removal of the packaging. Steriking® pouches are marked with a symbol indicating the correct peeling direction.

Where double packaging is required, it is important to position porous material against porous material and laminate against laminate because penetration of air and steam is only possible through the paper side. The inner packaging must not be folded so that the passage of sterilant remains unhindered.

### Limitations for Use

LT-Blueline Pouches and Rolls are not suitable for sterilization by hot, dry air or by steam.

### Sealing

All seals, including closure seal/seals should be smooth, i.e. without folds, bubbles, or wrinkles. The plastic film turns a darker shade where the sealing has taken place, making it easy to check that the seal is intact. Heat sealing devices must be capable of attaining the sealing conditions suitable

22.11.2022

Rev3.

for each specific sterile barrier system construction. Correct temperature, pressure and sealing time/speed combination need to be met. Preferably, only sealing devices manufactured and intended for medical use should be used.

The sealing temperature shall be in the range of 120°-130°C (248-266°F). As a recommendation the sealing force should be at minimum 95 N and at maximum 120 N. The sealing time of heat-sealer is not normally adjustable;



the speed of heat-sealer is usually from 8 m/min to 10 m/min. The needed parameters might vary due to the quality and type of the sealer used.

The seals need to be strong to withstand the sterilization process and handling yet providing a clean peel. Closing too strongly should be avoided. A manual test should be carried out for controlling the seal strength. Closures that compress the package should not be used, (e.g., ropes, strings, elastic bands, paperclips, staples or similar items).



### Labeling and writing

Writing or printing on pouches should only take place on the film side or on the Tyvek® outside the seals. Writing instrument should not have a potential for creating a hole or puncture in the sterile barrier system, i.e. ballpoint pens should not be used. Only markers intended for appropriate method of sterilization should be used. If labels are used they must not impede the sterilization process, i.e. must not block the breathable area of the package. Labels must not cover the seals or any necessary information such as LOT-number etc.



### Loading the sterilizer

Packages should be placed in the sterilizer loosely allowing the sterilizing agent to penetrate the pack. If a sterilization cycle must be repeated due to a malfunction or a cycle is aborted before completion, packages must be repacked prior to being placed into another sterilization cycle.

(HM)

The products are intended for single use only; product characteristics and performance cannot be guaranteed if sterilized more than once.

## STERRAD® Cycles

Compatible with following sterilization cycles:

NX Standard  
NX Advanced  
100NX Duo  
100NX Flex  
100NX Standard

## Device lumen dimensions:

Max lumen of 1 and 1 mm min ID for all cycles.

## Max Weights:

NX Standard Cycle:

75 mm x 200 mm, 100 mm x 250 mm, 150 mm x 300 mm, 160 mm x 440 mm, 160 mm x 600 mm:  
.088 pounds

205 mm x 390 mm, 250 mm x 390 mm, 250 mm x 500 mm:  
1.382 pounds

NX Advanced Cycle:

75 mm x 200 mm, 100 mm x 250 mm, 150 mm x 300 mm, 160 mm x 440 mm, 160 mm x 600 mm:  
.084 pounds

205 mm x 390 mm, 250 mm x 390 mm, 250 mm x 500 mm:  
1.393 pounds

100NX Duo Cycle:

75 mm x 200 mm, 100 mm x 250 mm, 150 mm x 300 mm, 160 mm x 440 mm, 160 mm x 600 mm:  
.104 pounds

205 mm x 390 mm, 250 mm x 390 mm, 250 mm x 500 mm:  
1.400 pounds

100NX Flex Cycle:

75 mm x 200 mm, 100 mm x 250 mm, 150 mm x 300 mm, 160 mm x 440 mm, 160 mm x 600 mm:  
.158 pounds

205 mm x 390 mm, 250 mm x 390 mm, 250 mm x 500 mm:  
1.400 pounds

100NX Standard Cycle:

75 mm x 200 mm, 100 mm x 250 mm, 150 mm x 300 mm, 160 mm x 440 mm, 160 mm x 600 mm: .158 pounds

205 mm x 390 mm, 250 mm x 390 mm, 250 mm x 500 mm:  
4.476 pounds

## Inspection

After sterilization each product is checked as to whether the packaging is intact, the changes communicated by indicator have taken place, and the product is clean and dry. In any unclear or ambiguous situation (e.g. wetness in pack) items should be considered as unsterile.

## Storage, Transportation and Shelf Life

It is recommended that products shall be kept in their original closed transportation carton and are stored in dry and clean conditions protected from direct sunlight and excessive moisture before they are taken into use.

It is recommended that the products are put to their end use within 5 years of manufacture. The recommended "Best before" date and the manufacturing date are stated on the carton label.

After sterilization the products are sorted for storage or delivery to the wards. The products should be stored in a dust-free place protected from sunlight, preferably in closed cabinets. It is recommended that the room climate has a humidity of 30 to 60% and a temperature of 15-25°C.



Maintenance of pack sterility is not only dependent on the packaging material and the method of sterilization: it is also dependent on handling, transport and storage conditions. The device is intended to allow sterilization of enclosed devices and to maintain sterility of the enclosed devices until used up to 3 years post sterilization. Any unnecessary handling of the packages should be avoided, as this would increase the risk of contamination. The level of protection is considerably enhanced by using a minimum of two layers - in other words by double pouching.

In cases where the transport or storage circumstances are particularly challenging, protective packaging such as a pouch made of impermeable multi-layer film can be used to protect the sterilized packages.

## Fiber-free opening and aseptic presentation

Follow the opening direction indicated on the pack. The seals on the upper corners of the packs should be unattached first. The package should then be opened by pulling the laminate away from the Tyvek® material slowly and evenly. Open only 1/3 - 2/3 of the packaging length. When opening large and / or heavy packages they need to be supported by a table or a tray.



Assistance may be needed to prevent any contamination of the packed instrument / instrument sets by accidentally touching the non-sterile outer surface of the packaging material. Double packaging ensures safe and

sterile opening. The inner pack remains sterile even on its outside until it is removed.

## Steriking LTS Flat Pouches

Art. Code	Size
LTS7520	75 x 200
LTS1025	100 x 250
LTS1530	150 x 300
LTS1644	160 x 440
LTS1660	160 x 600
LTS2038	205 x 390
LTS2538	250 x 390
LTS2550	250 x 500

## Serious Incidents

If any serious incident occurs in relation to the pouches or rolls, those should be reported to Wipak and the national competent authority in which user is established.

## Waste management

After use Steriking® sterilization pouches and rolls can be incinerated without producing toxic emissions. Of course, any contaminated product must be eliminated using a specialized method.



Steriking® is a registered trademark of Wipak.

Tyvek® is a registered trademark of DuPont.



*We care that you pack safely!*

# WIPAK

steriking@wipak.com · www.steriking.com