



healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
matt@healthmark.info

New Dry DuraSponge for Surgical Instruments

Fraser, MI, Aug 6, 2012- Healthmark Industries announces a new product for the cleaning of instrument that require immediate removal of soils after procedures, The Dry DuraSponge.

The Dry DuraSponge is the first step in the safe cleaning of sensitive instruments and fixtures, as well as all general surfaces. They retain 15% more of cleansing solution than standard sponges or paper wipes. The innovative tubular sponge is specifically designed for cleaning of flexible and rigid scopes.

Many companies offer sponges that are presoaked with enzyme detergents. They are often more expensive and usually not necessary since the user will combine them in an enzyme bath/basin which should contain the detergent mixed at the correct concentration. In some other cases, detergents aren't needed at all. Often flexible endoscopes only require water for patient bedside cleaning procedures. Check the device manufacturer's IFU before using.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.