



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
matt@healthmark.info

Healthmark Adds The Osteotome Rack to Product Line

Fraser, MI, Oct, 16, 2012 - Healthmark Industries announces the new Osteotome Rack to its Protech instrument trays line.

This purpose of this product is to protect osteotome tips. The innovative design keeps tips from damage during transportation and reprocessing. These racks make it easy and safe for handling in and out of the Operating Room. The racks are PTFE padded for complete edge protection with a fully flanged frame for a lightweight but durable design. Angled slots retain a maximum of 12 instruments in place when the bars are down.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.