



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
matt@healthmark.info

Healthmark Helps You Ensure Clean with TOSI®

Fraser, MI, Oct 10, 2012- Healthmark offers an effective device to provide the Sterile Processing Professional with a consistent, repeatable, reliable method for evaluating the cleaning effectiveness of their automated instrument washer.

The TOSI® is dried blood soil on a stainless steel coupon. It is analogous to dried blood on a stainless steel instrument. We mounted the plate in a plastic holder with a graduated gap. This is just like the areas on a surgical instrument which are blocked from direct spray action, such as box locks.

Failure to properly clean surgical instruments means that bioburden is "cooked" onto the instruments during sterilization - trapping contaminants where they can easily cause infection. Proper cleaning is critical - yet until the TOSI there was not been a device that consistently, and independently provides an objective challenge to cleaning effectiveness.

The TOSI® aligns perfectly with FDA, AAMI, AORN and other regulatory recommendations for a surrogate testing device: utilize the same kind of soil on the same type of surface and with the same kind of physical challenge as the instrument itself when reprocessed.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.

