



healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
matt@healthmark.info

Keep Your Endoscope Valves Safe

Fraser, MI, Dec. 7, 2012 - Healthmark Industries announces a new product to help keep endoscope valves and other detachable parts with the proper endoscope throughout reprocessing.

The Valvesafe™ locks onto the endoscope without concealing any part of the scope and is able to hold a maximum of four valves at once. This unique design made of polypropylene has large openings that allow detergents and disinfectants to flow through easily and effectively. The distinctive container allows valves and detachable parts to be kept with the parent endoscope. This will form a unique set of equipment that allows adequate drying in a contaminant free environment, while safe and simple storage keeps appropriate parts ready at all times.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.