



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
matt@healthmark.info

Healthmark Offers New USB Microscope to Optical Inspection Product Line

Fraser, MI, March 6, 2013 - Healthmark offers a new USB Microscope (MICR0007) that delivers improved image quality, and features which enhance the ability to closely inspect surgical instruments. This product helps ensure quality control for the CSSD by close visual inspection of reprocessed instruments.

The USB microscope has a 2M pixels LENS and CMOS sensor along with a magnification ratio of 5x-270x. The device comes with an aluminum stand which makes inspection much easier by freeing up both hands. The 8mm diameter design allows inspection in some cannulated items, including orthopedic shavers. The camera function used with the supplied software, allows you to capture still and video images when inspecting the instruments for documentation purposes.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.