



## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Stop Repeat Offenders with the Washer Rack Label**

Fraser, MI January 9, 2014 - Healthmark has announced a new label to their product line that are used to identify and track washer racks through the washing cycles when being tested with the TOSI™ and other ProFormance™ cleaning verification products that they offer.

This Washer Rack Label operates like a license plate for your washer rack, each label has a unique human and machine readable alphanumeric number. It will help determine washer racks that are repeatedly not performing to cleaning standards so the issue can be properly identified. Manufactured with anodized aluminum and robust adhesive backing, they are designed to withstand repeated washing cycles. They come in the size 2” x 0.625”. For more information about the Washer Rack Label, go to <http://hmark.com/racklabels.php?pmc=PR-WashRkLab>

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.