



**healthmark**  
INDUSTRIES CO.  
health care products

800-521-6224  
www.hmark.com

## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Healthmark Introduces New Green Mesh Valve Bag for Improved Endoscope Reprocessing**

Fraser, MI February 18, 2014 - Healthmark Industries announces the introduction of the new Green Mesh Valve Bag. This tough, lightweight, plastic mesh bag helps ensure the traceability of critical endoscope components from reprocessing through use by keeping the valves and scope together. Simply place the valves in the bag and then tie to the scope, thus securing the valves to their related flexible endoscope.

By using this product the time spent hunting in bins or drawers for the right valves will be reduced. Once in the procedure room, remove the valves and place on the scope and dispose of this inexpensive, single-use bag. For enhanced security and traceability add one of Healthmark's tamper evident seals and HangTime™ Labels. For more information go to <http://hmark.com/valvebag.php?pmc=GMB-PR>

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.