



## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Go Green with Made In America LED Magnifiers**

Fraser, MI, May 13, 2015 - Healthmark Industries has announced the addition of two new optical inspection products, Made In America LED Magnifiers. The 6" round LED magnifier has a 43" reach capability that adjoins to a screw down base, equipped with a standard 4 diopter (2x) lens and produces the highest lumen output in a lighted magnifier.

The 7" x 5 1/4" rectangular LED magnifier has a 43" reach that connects to a table edge clamp. Similar to the round magnifier this model has decreased glare, heat and energy consumption with multi-angle LEDs that make shadows obsolete. Assembled with low voltage (15-24v) power supplies, the magnifier auto shuts off every 4 hours to curb energy loss.

The Made In America LED Magnifiers produce a light output 8-watt LED array that is comparable to a 150-watt halogen. The LED technology produces more light using less energy. Their attractive design and vibrant colors will look good in any department.

For more information go to  
<http://www.hmark.com/USAmagnifiers.php?pmc=MIAMag-PR>  
or call 1-800-521-6224.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.