



**healthmark**  
INDUSTRIES CO.  
health care products

800-521-6224  
www.hmark.com

## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **New Silicone Heat Glove Offered by Healthmark**

**Fraser, MI, November 19, 2015** - Healthmark Industries has announced the Silicone Heat Glove to their instrument care glove line. Made from latex-free silicone, the 16.9 inch and 13.7 inch one-size fits-all Silicone Heat Glove is a conventional option for healthcare professionals performing tasks such as unloading steam sterilizers, washer-disinfectors or dealing with surface temperatures ranging from -40 °F to +464 °F. The protective heat glove can be worn ambidextrously by inserting either hand entirely into the glove with minimal effort. They are available in a shade of light blue, the reusable Silicone Heat Glove is suitable for reprocessing up to 273 °F.

For more information call 800-521-6224 or visit [www.hmark.com](http://www.hmark.com).

### **About Healthmark Industries Co., Inc.**

*Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, support and serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information or call 800-521-6224.*