



healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Check Endoscopes for Safety with the Now! Test

Fraser, MI, June 23rd 2016- Flexible endoscopes are a challenge to clean and disinfect/sterilize. Modern machines such as AER's are effective tools for rendering scopes safe for use on the next patient, but how do you determine if the process was effective? By using the new NOW! Test.

This simple and rapid test checks for Gram negative bacteria in about 12 hours, helping to ensure that it is safe to use on the next patient. Utilizing a unique enzyme detection method, the easy to read fluorometer checks for the Gram negative bacterial growth (<10 cfu) by reading telltale fluorescence in the recaptured water. If the fluorometer reading is positive for Gram negative bacteria, reprocess the endoscope following manufacturer guidelines prior to use.

Visit www.hmark.com or call 800-521-6224 for more information..

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, support and serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information or call 800-521-6224.