



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Healthmark Announces New Hanging Tag

Fraser, MI March 14th, 2017 - Healthmark Industries announced that the new Hanging Tag addition to their labeling product line.

Manufactured to be used with container systems, the white colored tag is 3.5 x 1.9 x 0.78 inches. The Hanging Tag features an ideal surface for barcode labels. Simply slide the tamper evident seal through cross-slot and affix to the container. The Hanging Tags are single-use, non-adhesive and heat resistant up to a temperature of 275° F.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.