



## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Healthmark Offers the Magic Touch**

Fraser, MI March 21st, 2017 - Healthmark has announced the addition of the Magic Touch Magnifier to their Optical Inspection line.

Designed for detailed inspection, the 3.75 x 3.75 x 1.8 Inch Magic Touch has a magnification of 3x and weighs 0.45 lbs. Equipped with three LED lights, the Magic Touch magnifier is a touch activated unit that features an on/off switch for uninterrupted lighting and twist-action focus that allows the user to zoom-in for greater detail. The magnifier includes two CR2016 batteries, a zip case and a microfiber lens cloth to help protect the glass lens from dust, debris and scratching.

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.