



## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Automatic Sampling for Channel Check**

Fraser, MI, May 24th, 2017 - Healthmark has announced the ACF-001 to its Endoscopy product line.

The ACF-001 is designed to inject specific volumes of water through channels of flexible endoscopes and other lumened devices when testing with the ChannelCheck™. It will not only automate the injection of water for sampling, but allow sampling of the entire suction/biopsy circuit of flexible endoscopes - even when laying horizontally. The ACF-001 comes with one solution source tube, a large and small delivery tube, as well as the following selectable sampling solution volumes: 7mL, 10mL and 20mL. The compact unit is equipped with a two-position selector switch for choosing the correct delivery tube line, and features an easy to navigate digital display screen. It also includes mountable hardware, a 250 mL graduated cylinder and an endoscope channel block.

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.