



## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

## **Healthmark Offers New Solutions for Wrapped Trays**

Fraser, MI, Sept. 18, 2017-Healthmark Industries has announced they are now offering Gorilla Bags™ as a new solution to help protect wrapped trays.

Made from spun polypropylene, the Gorilla Bags™ have a high tensile strength and are used as an outer protective layer applied over the standard hospital barrier packaging. The physical properties of the material offer increased resistance to tearing and splitting during transportation and storage. Due to the Gorilla Bags™ high permeability, they are applied before sterilization as it does not impede the steam from entering the package. Simply place the wrapped barrier package inside the Gorilla Bag™, fold the lip three times, seal the bag close with autoclave tape, label the package appropriately for further processing.

Visit [www.hmark.com](http://www.hmark.com) or call 800-521-6224 for more information.

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.