



## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

## **Healthmark Offers New Dry-Base® Sterisystem™ Trays**

Fraser, MI, Sept.28, 2017- Healthmark has announced the addition of Dry-Base® Sterisystem™ Trays and Lids to its ProTech® Instrument Tray product line.

Manufactured from highly polished AISI 304 grade stainless steel, the Dry-Base® Sterisystem™ Trays are designed to hold items during various stages of reprocessing and come outfitted with a flat base that remains in contact with condensation on wrap for better evaporation of condensate, as well as perforated side walls to facilitate effective processing of items. The stackable tray lids are equipped with four silicone corners, which are offered in four different colors: yellow, green, blue, and red. These silicone corners easily snap around the top wire of the tray to prevent small, lightweight items from falling out of the tray during processing. The Dry-Base® Sterisystem™ Trays and Lids can be paired with the Secur-Its™ for Wire Bottom Baskets from Healthmark Industries that simply attach and screw into the bottom of the tray for enhanced security of items during transportation and storage.

Visit [www.hmark.com](http://www.hmark.com) or call 800-521-6224 for more information.

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.