



FOR IMMEDIATE RELEASE:

New FIS-005 - The Next Generation of Flexible Inspection Scopes

Fraser, MI, May 1, 2018 - Healthmark has announced it launch of the FIS-005SK, the next generation of flexible inspections scopes to its ProSys™ Optical Inspection product line. Delivering unparalleled quality of image in an incredibly convenient and robust design. The scope is designed with a blue handle that includes a distal tip composed of a light source and camera lens at the end of a 110cm flexible blue shaft, which features white graduation marks. Manufactured for instruments 2.0mm in diameter or larger, the FIS-005 is the perfect tool to get a visualization of any potentially soiled or damaged item. With the Starter Kit Software is included, which installs on Windows 7, 8 & 10 PC's which powers the FIS-005 and allows for viewing, recording and documentation. The starter kit includes the scope, USB Drive and case.

Visit www.hmark.com or call 800-521-6224 for more information.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.