

FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith Healthmark Industries (586) 774-7600 msmith@hmark.com

New Flexible Inspection Scope That Enhances Visual Instrument Inspection

October 13, 2014 - Healthmark has announced the addition of the Flexible Inspection Scope to their Prosys Optical Inspection line of products.

The Flexible Inspection Scope includes a distal tip composed of a light source and camera lens at the end of a 50cm, flexible shaft. Designed for instruments 3.2mm in diameter or larger. The camera and light are powered by the USB connection on a PC.

Compatible with both Windows XP and Windows 7, the included software allows viewing and recording from most computers. Paired with the optional Flex Arm, the Flexible Inspection Scope can be securely fastened to a work stations to free both hands for manipulation of the scope and the target medical device. It is the perfect tool to visually inspect any device after cleaning, particularly those with internal channels and lumens.

For more information go to http://www.hmark.com/FIS.php?pmc=FIS-PR or call 1-800-521-6224.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.