



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Healthmark Announces New Window Style Tip Protector

May 16th, 2017 - Healthmark Industries has announced the new Window Style Tip Protector to their ProSys™ Instrument Care product line. Certain medical devices and instruments are a bit more delicate than others, and thus require equally delicate care and reprocessing practices.

The single-use Window Style Tip Protector is designed for protection of small sharp and fragile surgical instruments during the sterilization process. Manufactured from paper, the 50 x 125 mm Window Style Tip Protectors come with a 32mm wide slot for placement of instruments and features a transparent polyethylene pouch that secures the tip of the instrument in place to prevent from sliding around or out. The Window Style Tip Protectors are heat resistant up to 275° F and available for purchase in a box of 100.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.