


Brand Name of Product	Econo-Loks
Generic Name of Product	Tamper Evident Seals
Product Code Number(s)	7050 GN, 7050 RD
Intended Use:	To display evidence of tampering or use for access of critical equipment and supplies.
Range of Applications for Product	Crash carts, emergency boxes, cabinets, and other areas where demonstrating tampering or use is critical.
Key Specifications of Product	<ul style="list-style-type: none"> • Numbered • Sealed and broken by hand • Cannot be resealed once broken. • Variety of colors available.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	One hundred (100) Econo-Loks per box.
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Seal crash carts, emergency boxes, and other areas where evidence of tampering and use is critical.
Preparation for Use	N/A
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Pass arrow shackle through the opening of item to be sealed. 2. Firmly press arrow shackle into the locking device until it engages locking teeth. 3. Pull on seal to ensure it has been properly locked.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	Follow the policy of your facility for recording date, staff, and other information critical to management of critical inventory.
Special Warnings and Cautions	If a seal is broken, assume items have been tampered with, used, or removed and follow your facility's instructions for inspection and replacement of critical inventory.
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A

Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	Compatible with steam sterilization. Standard cycle: <ul style="list-style-type: none"> • 132 °C (270 °F) for 4 minutes • 135 °C (275 °F) for 3 minutes.
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	Security Products Brochure, Security Products Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com