

# Leak Tester Tester (LTT) Log Sheet

Department Name		Leak Tester Brand/Model		Serial Number	
-----------------	--	-------------------------	--	---------------	--

Date	Time	Staff Initials	LTT Results	Action if Pressure Out of Range

**Olympus MU-1 Pump:** >4.0 PSI

**Olympus MB-155 Cable:** Must be between 2.9-4.1 PSI

**Pentax Handheld/Connection:** >5.0 PSI

**Karl Storz Handheld/Connection:** >3.0 PSI

**Olympus handheld/connection:**

Unit of Pressure	Pressurize to	Correlation to LTT Gauge
kPa	25 units	3.6 psi on test fixture
mmHg	160 units	3.0 psi on test fixture
lb/in <sup>2</sup> (psi)	3 psi	3.0 psi on test fixture