

Instrument Trays

Robotic Instrument Cases

Scope Trays

Specialty Trays

Wire Baskets

Instru-Trays

Silicone Mats



Intelligent Solutions
for Instrument Care &
Infection Control


hmark.com

-  Cleaning Verification
-  Cleaning Chemistries
-  Protective Equipment
-  Labeling & Signage
-  Instrument Care
-  Sterilization Products
-  Storage Products
-  **Instrument Trays**
-  Security Products

da Vinci Arm Sterilization Case

Secure up to 8 robotic instruments


Manufactured of tough but lightweight anodized aluminum, this tray is configured to accommodate and secure robotic instrumentation for safe steam sterilization, transport and storage. The Secur-Its™ safely hold and protect up to 8 instruments with two 10 inch, 15 position universal brackets.



Robotic Arm Case

Secure 1 or 2 robotic instruments


Manufactured of tough but lightweight anodized aluminum, this tray is configured to accommodate and secure robotic instrumentation for safe steam sterilization, transport and storage. The Secur-Its™ safely hold and protect one or two instruments.



da Vinci Lens Case & Lid

Secure 1 or 2 robotic lenses

The da Vinci® Lens Case & Lid is manufactured of tough but lightweight anodized aluminum. Its designed for safe transportation and reprocessing of the da Vinci® instrument.




	Item Number	Description	Length	Width	Height	Din Size	Qty
A	31127 DV	da Vinci Arm Sterilization Case	27.00"	11.00"	3.50"	Extra Long	EA
B	300605-DVA	Robotic Arm Case	30.00"	6.00"	2.00"		EA
C	300605-DVC	da Vinci Lens Case & Lid	30.00"	6.00"	5.00"		EA

Bariatric Laparoscopic Tray

Two-level design holds up to 15 instruments

Large enough for the very long instruments used in bariatric surgeries, this double-decker system holds twice the number of instruments in the same amount of space. The laparoscopic insert includes graduated Secur-Its™ designed to accommodate up to 15 lap instruments.



Cysto Double-Decker Tray

Secure twice the number of instruments


The included Secur-Its™ can be placed virtually anywhere in this tray for maximum efficiency. The double-decker design holds twice the number of instruments in the same amount of space and the graduated Secur-Its™ are designed to hold instruments ranging from to 3 to 10mm in diameter.



Stainless-Steel Scope Trays

Safely secure delicate scopes


These tough but lightweight scope trays are an excellent alternative to factory supplied trays. Available in two sizes, each includes Secur-Its™ that can be mounted virtually anywhere they are needed within the tray.




	Item Number	Description	Length	Width	Height	Qty
A	51527	Bariatric Laparoscopic Tray	27.00"	11.00"	8.00"	EA
B	CYSTO-2214B	Cysto Double-Decker Tray	20.00"	10.00"	3.50"	EA
C	30111-SSC	Single Scope Tray	11.00"	3.00"	2.50"	EA
D	30118-SSC	Single Scope Tray	18.00"	3.00"	2.50"	EA

Clear-Top Scope Trays

Designed to secure rigid scopes



A



B

Clear-Top Scope Trays feature a lightweight, tough and durable construction. Available in two sizes, each option includes plug-in Secur-Its that are designed to snugly secure and protect two rigid scopes.

Long Scope Instrument Tray

Designed to secure long rigid scopes




C


Designed for long rigid scopes and other delicate instruments, this tray is manufactured of lightweight, durable polypropylene and is compatible with sterilization. Silicone grips securely hold the scope and instruments in place.

Bariatric Scope Trays

Designed to secure long bariatric scopes



D




E

The Bariatric Scope Trays are specially designed to hold the very long bariatric scopes. Healthmark offers a single-scope style and a multi-scope style. The multi-scope style is capable of storing 1-2 scopes and several instruments.


Customizable Full Size Trays

Multiple sizes for a variety of uses


Design a tray to meet your exact needs. Select from several sizes and customize with universal instrument brackets and silicone finger mats. Trays are made of anodized aluminum with perforations to allow for proper penetration of sterilant.




F



G



H



I

	Item Number	Details	Length	Width	Height	Din Size	Qty
	A	261615 Clear Top Scope Tray	16.00"	2.50"	1.50"		EA
	B	302215 Scope Tray	22.00"	3.00"	1.50"		EA
	C	SMT5510 Long Scope Instrument Tray	26.00"	4.25"	2.50"		EA
	D	30126 Bariatric Scope Tray - Single	26.00"	3.50"	3.50"		EA
	E	30630 Bariatric Scope Tray - Multi	30.00"	6.00"	2.50"		EA
	F	31421-H Universal Tray with Lid	21.00"	10.50"	2.50"	Full Din	EA
		31123-T Universal Tray with Insert & Lid	23.00"	11.00"	3.50"	Long	EA
	G	31123-T1 Universal Tray with Insert & Lid	23.00"	11.00"	5.00"	Long	EA
	H	31127-H Universal Tray with Lid	27.00"	11.00"	3.50"	Extra Long	EA
	I	30320H Universal Tray with Insert & Lid	20.00"	10.00"	3.50"	Full Din	EA



Metal Micro Protec Trays

Lightweight trays with smooth edges

These metal trays rival plastic trays in pricing, weight and smooth edges. Combine with silicone finger mats, or utilize the insert instrument holders to hold instruments in place. All are compatible with steam, EtO and Dry-Heat sterilization.



A



B




C



D




E



F



G



H



I



J

	Item Number	Details	Length	Width	Height	Qty
A	184014	Perforated Base and Cover	6.00"	2.50"	1.20"	EA
B	184016	Perforated Base and Cover	6.50"	4.00"	1.20"	EA
C	184019	Perforated Base and Cover	7.50"	4.00"	1.20"	EA
D	184021	Perforated Base and Cover	7.50"	4.00"	2.00"	EA
E	181453	Tray with Holder (Non-Perforated)	7.00"	11.00"	1.50"	EA
F	184453	Tray with Holder (Perforated)	7.00"	11.00"	1.50"	EA
G	181463	Base, Cover and Holder	7.00"	5.00"	1.25"	EA
H	182500	Stainless Steel Tray with 5 Clips	7.00"	3.33"	1.00"	EA
I	182062	Stainless Steel Cassette	7.00"	5.50"	1.00"	EA
J	182068	Twist-Locking Stainless Steel Cassette	11.00"	7.00"	1.00"	EA



Specialty Trays

Multiple sizes for a variety of uses


Plastic Instru-Trays are offered in a variety of sizes and colors. Each size is available with or without holes and they are ideal for various uses such as floor trays, mayo trays, custom trays, procedure trays or surgical trays.



A



B



C



D



E

	Item Number	Details	Length	Width	Height	Qty
A	30109HN	ER Hand Tray	10.00"	9.00"	3.00"	EA
B	30109CC	Camera Case with Brackets	10.00"	9.00"	3.00"	EA
C	30011KE	Key Elevator Box	10.00"	5.00"	2.00"	EA
	30011CU	Curretes Box	10.00"	5.00"	2.00"	EA
	30011RS	Rhoton Suction Box	10.00"	5.00"	2.00"	EA
	30011BP	Bipolar Set Box	10.00"	5.00"	2.00"	EA
D	SP-3014	Osteotome Tray	11.00"	9.00"	3.00"	EA
E	31020-H	Protech Base, Cover & Insert	20.00"	10.00"	2.625"	EA



Hinged-Top Micro Trays

Silicone Mats & Secur-Its sold separately

Designed to organize micro-surgical instruments & hardware, Hinged-Top Micro Trays include a hinged cover and snap latch to allow fast and secure operation. These trays are designed for use in all forms of sterilization.



Specialty Hinged-Top Micro Trays

For Scopes, Cameras and Arthroscopes

Specialty Hinged-Top Micro Trays offer the same features as the original Hinged-Top Micro Trays but they include Silicone Plug-In Secur-Its designed to help secure instruments such as scopes, cameras and arthroscopes.




	Item Number	Details	Length	Width	Height	Qty
A	30732-H	Hinged-Top Tray - Perforated	7.25"	3.25"	1.125"	EA
B	30842-H	Hinged-Top Tray - Perforated	8.75"	4.50"	1.75"	EA
C	31042-H	Hinged-Top Tray - Perforated	10.25"	4.50"	1.75"	EA
D	31055-H	Hinged-Top Tray - Perforated	10.00"	5.00"	3.25"	EA
E	31462-H	Hinged-Top Tray - Perforated	14.00"	6.25"	2.50"	EA
F	31772-H	Hinged-Top Tray - Perforated	17.50"	7.75"	2.375"	EA
G	30732	Hinged-Top Tray - Non-Perforated	7.25"	3.25"	1.125"	EA
H	30842	Hinged-Top Tray - Non-Perforated	8.75"	4.50"	1.75"	EA
I	31042	Hinged-Top Tray - Non-Perforated	10.25"	4.50"	1.75"	EA
J	31055	Hinged-Top Tray - Non-Perforated	10.00"	5.00"	3.25"	EA
K	31462	Hinged-Top Tray - Non-Perforated	14.00"	6.25"	2.50"	EA
L	31772	Hinged-Top Tray - Non-Perforated	17.50"	7.75"	2.375"	EA
M	50Z900	Hinged-Top Tray - Perforated	8.00"	1.75"	1.75"	EA
N	31042-22	Hinged-Top Arthro Tray	10.25"	4.50"	1.75"	EA
O	31772-22	Economy Scope Tray	17.50"	7.75"	2.375"	EA
P	31462-22	Camera Tray	14.00"	6.25"	2.50"	EA




Clear-Top Micro Protech Trays

Designed for microsurgical instruments


Organize microsurgical instruments and accessories in these customizable trays. Removable, clear, snap-on covers make for easy identification of contents. All trays are compatible with Steam, EtO, and low temperature sterilization.




R




S




T




U



A



B



C




D




E




F




G




H



I




J




K




L




M




N



O



P



Q

	Item Number	Details	Length	Width	Height	Qty	
	A	0278	Base & Cover	2.70"	1.50"	0.75"	EA
	B	02781	Base & Cover	2.70"	1.50"	1.00"	EA
	C	027814	Base & Cover	2.70"	1.50"	1.30"	EA
	D	0278-M	Silicone Mat for 0278 Series Trays	2.70"	1.50"		EA
	E	0625	Base & Cover	6.00"	2.50"	0.80"	EA
	F	062514	Base & Cover	6.00"	2.50"	1.30"	EA
	G	0625-M	Silicone Mat for 0625 Series Trays	6.00"	2.50"		EA
	H	6540	Base & Cover	6.50"	4.00"	0.75"	EA
	I	6540-M	Silicone Mat for 6540 Tray	6.50"	4.00"		EA
	J	100675	Base & Cover	10.00"	6.00"	0.75"	EA
	K	100675-M	Silicone Mat for 100675 Trays	10.00"	6.00"		EA
	L	7540	Base & Cover	7.50"	4.00"	0.75"	EA
	M	754015	Base & Cover	7.50"	4.00"	1.50"	EA
	N	7540-M	Silicone Mat for 7540 Series Trays	7.50"	4.00"		EA
	O	100615	Double-Decker Base & Cover	10.00"	6.00"	1.50"	EA
	P	100615B	Full Insert Tray for 100615				EA
	Q	100615-M	Silicone Mat for 100615 Trays	10.00"	6.00"		EA
	R	101525	Double-Decker Base & Cover	15.00"	10.00"	1.50"	EA
	S	101525B	Full Insert Tray for 101525	15.00"	10.00"		EA
	T	101525-M	Silicone Mat for 100615	15.00"	10.00"		EA
	U	101525B-M	Silicone Mat for 100615B	15.00"	10.00"		EA



Wire Baskets

Multiple sizes for a variety of uses


Wire baskets are designed to stack one atop the other for convenient and safe storage of instruments. Sizes are compatible with European-made closed container systems and washer-decontaminators.



A



B




C



D



E



G



H




I

Fine Mesh Baskets with Snap Closure

Secure small items in instrument washers

Manufactured from stainless steel, these small baskets are ideal for safely washing small items that often get lost in instrument washers. The small perforations ensure effective water penetration during reprocessing.



J



K

	Item Number	Description	Length	Width	Height	Din Size	Qty
A	75.0005.4	Wire Basket	9.50"	9.875"	2.00"	Half Din	EA
B	75.0005.5	Wire Basket	9.50"	9.875"	2.75"	Half Din	EA
C	75.0005.6	Wire Basket	18.875"	9.875"	4.00"	Full Din	EA
D	75.0005.9	Wire Basket	18.875"	9.875"	2.00"	Full Din	EA
E	75.0194.2	Cover for Wire Baskets	18.875"	9.875"			EA
F	75.0091.2	Cover for Wire Baskets	18.875"	9.875"			EA
G	75.0194.4	Wire Basket	18.875"	9.875"	2.75"		EA
H	72.1400.5	Washer Basket	23.64"	15.76"	3.94"		EA
I	72.1400.6	Washer Basket	23.64"	15.76"	7.09"		EA
J	CIL-004	Fine Mesh Basket	105mm	70mm	25mm		EA
K	CIL-005	Fine Mesh Basket	80mm	40mm	20mm		EA



Perforated Stainless Steel Trays

Durability to last a lifetime

The perforated design ensures maximum steam penetration but limits the chance for delicate instrument points to become damaged. They are made with a strong stainless steel frame and ensure a construction to last a lifetime.



Item Number	Description	Length	Width	Height	Din Size	Qty
A	PSST-001	15.50"	5.00"	3.50"	Three-Quarter Din Narrow	EA
B	PSST-002	20.00"	5.00"	3.50"	Full Din Narrow	EA
C	PSST-003	20.25"	13.00"	3.50"	Wide Body	EA
D	PSST-004	25.00"	15.00"	3.50"	Extra Long	EA
E	PSST-005	10.50"	15.00"	3.50"	Three-Quarter Din	EA
F	PSST-006	30.00"	6.00"	3.50"	Extra Long Narrow	EA
G	PSST-007	10.50"	10.00"	3.50"	Half Din	EA
H	PSST-008	26.00"	9.00"	3.50"	Extra Long	EA
I	PSST-009	10.00"	6.50"	3.50"	Quarter Din	EA
J	PSST-010	10.50"	20.00"	3.50"	Wide Body	EA
K	PSST-011	24.00"	6.00"	3.50"	Extra Long Narrow	EA
L	PSST-012	20.25"	13.00"	3.50"	Wide Body	EA
	PSST-013	25.25"	8.00"	3.50"	Extra Long	EA



High-Sided Trays

Multiple sizes for a variety of uses

High-Sided Trays are made of durable, non-corrosive surgical stainless steel. Available perforated and non-perforated, they are designed to withstand the demanding environment of surgery and surgical instrument reprocessing.




Item Number	Description	Length	Width	Height	Qty
A	HDO-0952	8.875"	5.00"	2.00"	EA
B	HDO-0952-NP	8.875"	5.00"	2.00"	EA
C	HDO-1062	10.00"	6.50"	2.00"	EA
D	HDO-1062-NP	10.00"	6.50"	2.00"	EA
E	HDO-1210	12.75"	10.375"	2.50"	EA
F	HDO-1272	12.25"	7.75"	2.25"	EA
G	HDO-1610	16.50"	10.00"	2.50"	EA
H	HDO-1610-NP	16.50"	10.00"	2.50"	EA
I	HDO-2012	20.75"	12.75"	2.50"	EA



Stainless Steel Instru-Trays

Multiple sizes for a variety of uses

Stainless Steel Instru-Trays are offered in a variety of sizes. Each size is available with or without holes and they are ideal for various uses such as floor trays, mayo trays, custom trays, procedure trays or surgical trays.



	Item Number	Style	Material	Length	Width	Height	Qty
A	HCO-1065	Perforated	Stainless Steel	10.00"	6.50"	0.75"	EA
B	HCO-1309	Perforated	Stainless Steel	13.688"	9.813"	0.75"	EA
C	HCO-1309-NP	Non-Perforated	Stainless Steel	13.688"	9.813"	0.75"	EA
D	HCO-1510	Perforated	Stainless Steel	15.125"	10.625"	0.75"	EA
E	HCO-1912	Perforated	Stainless Steel	19.125"	11.625"	0.75"	EA
F	HCO-1912-NP	Non-Perforated	Stainless Steel	19.125"	11.625"	0.75"	EA
G	HCO-1711	Perforated	Stainless Steel	17.00"	11.625"	0.75"	EA
H	HCO-2116	Perforated	Stainless Steel	21.25"	16.25"	0.75"	EA



Plastic Sterilization Baskets

Multiple sizes & colors for a variety of uses

Unlike metal, these durable polypropylene plastic trays have no sharp edges to tear wrap and they are easy to handle and transport. They are ideal for the sterilization of all hospital instruments and supplies and are autoclavable to 285°F.



	Item Number	Description	Length	Width	Height	Din Size	Qty
A	PIT2015	Plastic Sterilization Basket	8.00"	6.00"	2.00"		EA
B	DINT2525	Plastic Sterilization Basket	10.00"	10.00"	2.00"		EA
C	PIT2715	Plastic Sterilization Basket	10.64"	5.91"	1.18"		EA
D	MIT30251	Plastic Sterilization Basket	11.80"	9.875"	2.00"		EA
E	MIT40301	Plastic Sterilization Basket	16.50"	12.00"	3.00"		EA
F	DINT4825	Plastic Sterilization Basket	20.00"	10.00"	2.00"		EA
G	1511-7	Plastic Sterilization Peg Rack	19.75"	19.75"	3.75"		EA
H	DT2-1218	Plastic Sterilization Basket	16.00"	10.00"	3.00"	Three-Quarter Din	EA
I	DT2-2136	Plastic Sterilization Basket	17.75"	10.75"	4.00"		EA
J	1511-4	Plastic Sterilization Basket	15.125"	11.25"	4.375"		EA


Plastic Instru-Trays

Multiple sizes & colors for a variety of uses

Plastic Instru-Trays are offered in a variety of sizes and colors. Each size is available with or without holes and they are ideal for various uses such as floor trays, mayo trays, custom trays, procedure trays or surgical trays.




A




B




C



D



E



F



G



H




I



K



J



L


	Item Number	Style	Color(s)	Length	Width	Qty	
	A	41520	Non-Perforated	Blue	20.00"	15.00"	EA
	B	41520H	Perforated	Blue	20.00"	15.00"	EA
	C	41317	Non-Perforated	Blue, Gray, Red	17.75"	13.813"	EA
	D	41317H	Perforated	Blue, Gray, Red	17.75"	13.813"	EA
	E	41116	Non-Perforated	Blue, Gray, Mauve, Red	16.00"	11.75"	EA
	F	41116H	Perforated	Blue, Gray, Mauve, Red	16.00"	11.75"	EA
	G	41013	Non-Perforated	Blue, Gray, Mauve, Red	13.563"	10.438"	EA
	H	41013H	Perforated	Blue, Gray, Mauve, Red	13.563"	10.438"	EA
	I	4106	Non-Perforated	Blue, Green, Pink, Purple	10.125"	6.625"	EA
	J	4106H	Perforated	Blue, Green, Pink, Purple	10.125"	6.625"	EA
	K	4116	Non-Perforated	Blue	9.00"	6.00"	EA
	L	4116H	Perforated	Blue	9.00"	6.00"	EA



Silicone Finger Mats

Protect delicate surgical instruments

Protect delicate surgical instruments with these robust silicone mats. Silicone Finger Mats gently but firmly cushion instruments from damage and help to hold them in place, preventing damage. Choose from a variety of sizes.




Item Number	Description	Color	Length	Width	Qty
2-2025	Silicone Finger Mat	Blue	7.75"	7.75"	EA
2-1405	Silicone Finger Mat	Blue	7.80"	3.40"	EA
A 2-3215	Silicone Finger Mat	Blue	13.25"	8.25"	EA
2-1055	Silicone Finger Mat	Blue	17.50"	7.50"	EA
2-1215	Silicone Finger Mat	Blue	18.30"	8.35"	EA
2-4005	Silicone Finger Mat	Blue	19.75"	9.40"	EA
2-6005	Silicone Finger Mat	Blue	21.10"	9.80"	EA
6-1914	Silicone Finger Mat	Blue	19.00"	14.00"	EA
2111-M	Silicone Finger Mat	Blue	21.00"	11.00"	EA
307-M	Silicone Finger Mat	Blue	7.25"	3.25"	EA
308-M	Silicone Finger Mat	Blue	8.75"	4.5"	EA
3100-M	Silicone Finger Mat	Blue	10.00"	10.00"	EA
310-M	Silicone Finger Mat	Blue	10.00"	4.00"	EA
314-M	Silicone Finger Mat	Blue	14.00"	6.25"	EA
317-M	Silicone Finger Mat	Blue	17.50"	7.75"	EA
SP2313	Silicone Finger Mat	Pink	9.00"	5.00"	EA
SP4622	Silicone Finger Mat	Pink	9.00"	8.75"	EA
SP2719	Silicone Finger Mat	Pink	10.875"	7.625"	EA
B SR10	Silicone Finger Mat	Pink	39.40"	9.46"	EA



Brackets for Customizable Trays

Easily apply to a variety of trays

Add brackets to Healthmark's wide variety of customizable trays. Select the bracket styles that work best for the instruments to be secured, and easily apply them to a ProTech™ tray with the provided screw for a custom solution.



Secur-Its for Wire Bottom Baskets

Prevent instruments from falling out of place

These Secur-Its™ safe-guard your medical instruments to prevent them from falling out of place and they are heat-resistant up to 275 °F. Simply arrange and screw into the bottom of the tray to accommodate the instrument size.



Item Number	Description	Qty
A P-01	Post	EA
B TRI-CUT BRCKT	Universal Tri-Cut Bracket (3, 5 & 10mm) - 1.00"	EA
C 56-UB	6 Position Universal Bracket - 5.00"	EA
D 710-UB	10 Position Universal Bracket (3, 5 & 10mm) - 7.50"	EA
E 1015-UB	15 Position Universal Bracket - 10.00"	EA
F 7501542	Secur-Its for Wire Bottom Baskets - 12mm Dia.	PKG-4
G 7501544	Secur-Its for Wire Bottom Baskets - 12.5mm, 6.7mm & 4.3mm Dia.	PKG-4
H 7501545	Secur-Its for Wire Bottom Baskets - 5mm, 3mm & 1.8mm Dia.	PKG-4



Silicone Plug-In Secur-Its™

Easily apply to a variety of trays

Add Silicone Plug-In Secur-Its™ to Healthmark's wide variety of customizable trays. Select the bracket styles that work best for the instruments to be secured, and easily apply them to a ProTech™ tray by simply popping them in.



A




B




C



D



E



F



G



H



I



J



K



L



M



N



O



P




Q



Magnetic Mat

Prevent items from falling from surfaces

Designed with strong round-raised magnets, this mat is intended to be used as a hands-free "neutral zone" sharps transfer area. This reusable, waterproof mat is heat-resistant up to 275°F, conforms to all surface shapes and can be fastened to drapes with a towel clamp. It also includes cut-lines that can be split in half to create two smaller sizes.



R

	Item Number	Description	Slot Size(s)	Number of Slots	Number of Plugs	Qty	
	A	A237591K	Plug-In Secur-It	4mm	1	2	EA
	B	A237592K	Plug-In Secur-It	4mm	2	2	EA
	C	A237593K	Plug-In Secur-It	4mm	3	2	EA
	D	A237593	Plug-In Secur-It	1mm	3	2	EA
	E	A237596	Plug-In Secur-It	1mm	6	2	EA
	F	A387596	Large Plug-In Secur-It	1mm	6	3	EA
	G	A3875912	Large Plug-In Secur-It	1mm	12	3	EA
	H	A387592	Large Plug-In Secur-It	4mm	2	3	EA
	I	A387593	Large Plug-In Secur-It	4mm	3	3	EA
	J	A387594	Plug-In Secur-It	4mm	4	3	EA
	K	A387595	Plug-In Secur-It	4mm	5	3	EA
	L	A2616901	Plug-In Secur-It	12mm, 8mm	2	2	EA
	M	A2616902	Plug-In Secur-It	8mm	2	2	EA
	N	A32759P2	Plug-In Secur-It	5mm	1	3	EA
	O	A327593	Plug-In Secur-It	5mm	3	3	EA
	P	A327593L	Plug-In Secur-It	12mm, 5mm	3	3	EA
	Q	A32759P1	Plug-In Secur-It	12mm	1	3	EA
	R	11116	Magnetic Mat, 16" x 12"				PK-2



A Getinge company

Intelligent Solutions
for Instrument Care &
Infection Control

hmark.com

Healthmark, A Getinge company

18600 Malyn Blvd.

Fraser, Michigan 48026

USA

Call Us

Toll Free: 800.521.6224

International: +1.586.774.7600

Fax: 586.491.2113

Email Us

Contact: healthmark@hmark.com

Order Inquiries: hmorders@hmark.com

