

Storage Products

Tote Boxes

Open-Front Containers

Utility Carts

Storage Accessories









Disposable Cart Covers

Reusable Cart Covers



Intelligent Solutions
for Instrument Care &
Infection Control

hmark.com



-  Cleaning Verification
-  Cleaning Chemistries
-  Protective Equipment
-  Labeling & Signage
-  Instrument Care
-  Sterilization Products
-  Storage Products
-  Instrument Trays
-  Security Products



Open-Front Containers

Stackable bins in a variety of sizes & colors

Contents are easily accessible and easily identified through large open fronts. Space can be maximized by stacking bins one-on-top the other. Slide your own labels into large label slots. Accessories include dividers and covers for select sizes.



	Item Number	Color	Length	Width	Height	Divider	Cover	Qty
A	30-210 BL	Blue	5.375"	4.125"	3.000"	40-210		EA
B	30-210 RD	Red	5.375"	4.125"	3.000"	40-210		EA
C	30-210 YL	Yellow	5.375"	4.125"	3.000"	40-210		EA
D	30-220 BL	Blue	7.375"	4.125"	3.000"	40-220	30-221	EA
E	30-220 RD	Red	7.375"	4.125"	3.000"	40-220	30-221	EA
F	30-220 YL	Yellow	7.375"	4.125"	3.000"	40-220	30-221	EA
G	30-220 GREEN	Green	7.375"	4.125"	3.000"	40-220	30-221	EA
H	30-239 BL	Blue	10.750"	8.250"	7.000"	40-239		EA
I	30-239 RD	Red	10.750"	8.250"	7.000"	40-239		EA
J	30-239 YL	Yellow	10.750"	8.250"	7.000"	40-239		EA
K	30-239 GREEN	Green	10.750"	8.250"	7.000"	40-239		EA
L	30-230 BL	Blue	10.750"	11.000"	5.000"	40-230		EA
M	30-230 RD	Red	10.750"	11.000"	5.000"	40-230		EA
N	30-230 YL	Yellow	10.750"	11.000"	5.000"	40-230		EA
O	30-230 GREEN	Green	10.750"	11.000"	5.000"	40-230		EA
P	30-240 BL	Blue	10.875"	5.500"	5.000"	40-245		EA
Q	30-240 RD	Red	10.875"	5.500"	5.000"	40-245		EA
R	30-240 YL	Yellow	10.875"	5.500"	5.000"	40-245		EA
S	30-240 GREEN	Green	10.875"	5.500"	5.000"	40-245		EA
T	30-235 BL	Blue	14.750"	8.250"	7.000"	40-230		EA
U	30-235 RD	Red	14.750"	8.250"	7.000"	40-230		EA
V	30-235 YL	Yellow	14.750"	8.250"	7.000"	40-230		EA
W	30-235 GREEN	Green	14.750"	8.250"	7.000"	40-230		EA



Open-Front Containers

Stackable bins in a variety of sizes & colors

Contents are easily accessible and easily identified through large open fronts. Space can be maximized by stacking bins one-on-top the other. Slide your own labels into large label slots. Accessories include dividers and covers for select sizes.



A



B



C



D



E



F



G



H



I



J



K



L



M



N



O



P



Q



R



S

	Item Number	Color	Length	Width	Height	Divider	Qty	
	A	30-250 BL	Blue	14.75"	16.50"	7.00"	40-245	EA
	B	30-250 RED	Red	14.75"	16.50"	7.00"	40-245	EA
	C	30-250 YL	Yellow	14.75"	16.50"	7.00"	40-245	EA
	D	30-250 GREEN	Green	14.75"	16.50"	7.00"	40-245	EA
	E	30-265 BL	Blue	18.00"	8.25"	9.00"	40-265	EA
	F	30-265 RD	Red	18.00"	8.25"	9.00"	40-265	EA
	G	30-265 YL	Yellow	18.00"	8.25"	9.00"	40-265	EA
	H	30-270 BL	Blue	18.00"	16.50"	11.00"	40-270	EA
	I	30-270 RD	Red	18.00"	16.50"	11.00"	40-270	EA
	J	30-270 YL	Yellow	18.00"	16.50"	11.00"	40-270	EA
	K	30-270 GREEN	Green	18.00"	16.50"	11.00"	40-270	EA
	L	900-00012	Blue	23.875"	16.50"	11.00"		EA
	M	QUS260	Blue	18.00"	11.00"	10.00"		EA
	N	900-00013	Blue	23.875"	8.25"	9.00"		EA
	O	900-00014	Blue	11.00"	5.50"	5.00"		EA
	P	QSB206	Blue	23.625"	6.625"	6.000"		EA
	Q	QSB104 BL	Blue	17.875"	6.625"	4.00"	DIV-104	EA
	R	900-00010	Blue	11.625"	4.125"	4.00"		10
	S	900-00011	Blue	11.625"	6.625"	4.00"		EA



Ribbed Tote Boxes

Dividable plastic containers

Regain space with compartmentalized ribbed containers and comply with The Joint Commision regulations by removing contents from outside shipping boxes and storing in impact-resistant containers that can be cleaned & disinfected.



A



B



C



D



E



F



G



H



I



J



K



L



	Item Number	Description	Length	Width	Height	Short Divider	Long Divider	Cover	Qty
A	MTB076	Ribbed Tote Box	10.750"	10.750"	6.125"			07C	EA
B	200-5 GY	Ribbed Tote Box	11.000"	8.250"	5.000"				EA
C	MTB086	Ribbed Tote Box	13.000"	13.000"	6.125"				EA
D	MTB096	Ribbed Tote Box	15.000"	15.000"	6.125"	DIV-096	DIV-096		EA
	300-2 GY	Ribbed Tote Box	16.375"	10.875"	2.500"	DIV-321	DIV-322		EA
E	300-4 GY	Ribbed Tote Box	16.375"	10.875"	4.250"				EA
	300-5 GY	Ribbed Tote Box	16.375"	10.875"	5.125"				EA
	300-6 GY	Ribbed Tote Box	16.375"	10.875"	6.000"				EA
	300-7 GY	Ribbed Tote Box	16.375"	10.875"	7.000"	DIV-361	DIV-362		EA
	300-8 GY	Ribbed Tote Box	16.375"	10.875"	8.000"				EA
	16-4 GY	Ribbed Tote Box	17.750"	5.500"	4.000"		DIV-5042		EA
	33-6 GY	Ribbed Tote Box	17.750"	10.875"	6.000"	DIV-3361			EA
F	50-3 GREY	Ribbed Tote Box	17.750"	16.375"	3.000"	DIV-5021	DIV-5022	50R CL	EA
	50-5 GY	Ribbed Tote Box	17.750"	16.375"	5.000"				EA
	50-7 GY	Ribbed Tote Box	17.750"	16.375"	7.000"				EA
G	MTB01804	Ribbed Tote Box	18.250"	15.125"	4.000"	DIV-1814	DIV-1824	1804C	EA
	MTB02106	Ribbed Tote Box	19.250"	13.500"	6.000"	DIV-2116			EA
H	MTB02109	Ribbed Tote Box	19.250"	13.750"	9.000"	DIV-2119			EA
I	MTB02006	Ribbed Tote Box	20.500"	14.500"	6.000"	DIV-2016	DIV-2026		EA
	MTB181604	Ribbed Tote Box	21.000"	19.000"	4.250"				EA
	MTB181608	Ribbed Tote Box	21.000"	19.000"	8.750"	DIV-1618			EA
J	400-4 GY	Ribbed Tote Box	22.375"	17.500"	4.250"				EA
K	400-5 GY	Ribbed Tote Box	22.375"	17.500"	5.000"				EA
	400-6 GY	Ribbed Tote Box	22.375"	17.500"	6.000"	DIV-451	DIV-452		EA
	400-7 GY	Ribbed Tote Box	22.375"	17.500"	7.000"				EA
	400-8 GY	Ribbed Tote Box	22.375"	17.500"	8.000"				EA
L	MTB01210	Ribbed Tote Box	27.000"	20.375"	9.000"				EA
	MTB02805	Ribbed Tote Box	29.250"	15.000"	4.500"				EA


Additional sizes available.



Ribless Tote Boxes

Several sizes for transportation or storage

Transport clean or dirty equipment and save space by using containers of the right size for efficient storage. Made of tough but lightweight plastics. Easy to clean with no nooks or crannies. Snap-on covers to keep supplies safe and dust-free.




	Item Number	Description	Length	Width	Height	Cover	Qty
A	GP110801	Ribless Tote Box	10.375"	7.875"	1.125"		EA
	GP100504	Ribless Tote Box	10.875"	6.000"	3.875"		EA
B	GP121004	Ribless Tote Box	12.125"	10.500"	4.000"	GP1210P	EA
C	GP151106	Ribless Tote Box	14.750"	11.500"	6.000"		EA
D	GP151505	Ribless Tote Box	14.750"	14.750"	5.000"		EA
	GP161104	Ribless Tote Box	15.750"	10.875"	3.875"		EA
E	GP181204	Ribless Tote Box	17.750"	11.750"	3.875"		EA
	GP181505	Ribless Tote Box	18.500"	15.000"	4.500"		EA
F	GP180804	Ribless Tote Box	18.875"	8.500"	4.000"		EA
G	GP201006	Ribless Tote Box	19.250"	9.250"	6.250"	GPL2010C	EA
H	GP191506	Ribless Tote Box	19.250"	15.000"	6.250"		EA
I	GP201605	Ribless Tote Box	20.250"	16.250"	5.500"	GPL2016C	EA
J	GP211304	Ribless Tote Box	20.750"	13.375"	3.875"		EA
	GP211705	Ribless Tote Box	21.000"	17.000"	5.000"		EA
K	GP210803	Ribless Tote Box	21.125"	8.125"	3.500"		EA
L	GP211406	Ribless Tote Box	21.500"	14.250"	6.000"		EA
M	GP221805	Ribless Tote Box	22.750"	18.625"	5.500"		EA
	GP230806	Ribless Tote Box	23.500"	8.125"	6.000"		EA
	GP251807	Ribless Tote Box	24.875"	18.500"	6.750"		EA
N	GP251910	Ribless Tote Box	25.500"	19.750"	10.000"		EA
O	GP261204	Ribless Tote Box	27.375"	12.000"	4.000"		EA

Flip-Top Containers

For safe storage & transport of supplies

Flip-Top Containers have attached, hinged covers for convenient and safe storage and transport of supplies. They provide compact storage by nesting one in another when not in use. The lockable tops offer security and use control. The molded in handle provides secure handling.



A

Replacement Patient Bins

Label holder helps identify contents

Replace old, discolored and damaged picking station bins with new ones. These tan Replacement Patient Bins are manufactured of tough polystyrene plastic and include a molded in label holder for identifying contents.



B

C

Open Front Wire Basket

Conveniently store supplies

Conveniently store supplies in Open Front Wire Baskets. A variety of supplies are easily identified with the see-through baskets. Open Front Wire Baskets can be easily hooked to a shelf, with no hardware required.



D

Infinite Divider Box

Dividers can be placed anywhere

Imagine the arrangements you can create with dividers that can be positioned virtually anywhere. Wasted space is eliminated because there are no fixed divider locations. These trays are made of super strong translucent polypropylene.




E

Stainless Steel Peg Boards

The holes accept a variety of accessories

The Stainless Steel Peg Boards are a great solution for wall-mounted storage. They include two 16" x 32" panels and include matching stainless screws. They work well for a variety of applications and help to induce a clean environment. They have holes that accept a wide range of hooks, bins, shelves and other accessories.



F

Undershelf Slides

Increase bin storage space on wire shelves

Undershelf Slides are stainless steel brackets that are designed to be mounted underneath wire shelves and hold compatible hanging plastic bins to add bin storage on wire shelving units.



G

	Item Number	Details	Length	Width	Height	Qty
	QDC2115-12	Flip-Top Container	21.50"	15.25"	12.75"	EA
A	QDC2420-12	Flip-Top Container	24.00"	20.00"	12.50"	EA
B	900-00015	Replacement Patient Bin	11.00"	2.75"	2.50"	EA
C	900-00016	Replacement Patient Bin	11.00"	2.75"	3.3125"	EA
D	75.0294.5	Open Front Wire Basket	13.75"	8.25"	7.875"	EA
E	T600	Infinite Divider Box	11.00"	6.625"	2.50"	EA
	T900	Infinite Divider Box	13.75"	8.75"	1.75"	EA
F	ALGBRD1632STLS	Stainless-Steel Peg Board	16.00"	32.00"		2
G	US-18	Undershelf Slide	18.00"			2
	US-21	Undershelf Slide	18.00"			2

EndoDolly

Designed to hold up to two scopes per pole

The EndoDolly™ is designed for hanging scopes and stands on a five-wheel locking base for ease of transportation. It is equipped with 3 self-adjusting extension poles that can be individually raised or lowered by simply touching a pressure button. The user can extend the 2 outside poles to a maximum height of 83.5 inches or lowered to a minimum height of 52.4 inches, with the center pole able to be adjusted to a maximum height of 107 inches or lowered to a minimum height of 58.7 inches. Located at the top of each pole is a lock tight hook holder that can accommodate two scopes per pole of any length and supports flexible scopes up to 15 lbs. per pole.



EndoWally

Designed to hang scopes with a wall mount


The EndoWally™ is designed for hanging scopes and is built with wall mount support arms that pivot to allow the user to swing the scope into place as needed. The EndoWally™ is equipped with a self-adjusting pole that can be raised or lowered by simply touching a pressure button. For stability, the user should ensure that 4 inches of the pole is always sticking out through the bottom support arm, which means the max the pole can be extended to is 61 inches. Located at the top of the pole is a lock tight hook holder that can accommodate two scopes of any length and has a weight capacity of 40 lbs.



Scope Holder

Designed to hang scopes vertically

Crafted from Plexiglas®, the Scope Holder is a wall-mountable design manufactured for assistance in testing various scopes prior to medical procedures. The vertical storage unit provides medical professionals with additional space, while keeping the scope in the proper upright position for thorough analysis before use.



Louvered Panel

Save space with wall-mounted storage

The Louvered Panel is a great solution for wall-mounted storage. For use with Open Front Containers of various sizes and colors, allowing the user to create a custom organization system. The Louvered Panel is made of steel and is powder coated gray. The load capacity is 175 lbs.



Wall-Mounted Bin Rail

Mount a variety of open front containers

This 48" long x 3" high steel rail can be mounted to any surface or workbench and is designed to hold multiple Open Front Containers in a variety of sizes. The Wall-Mounted Bin Rail has a load capacity of 125 lbs.



Item Number	Description	Length	Width	Qty
A	110405-SP EndoDolly			EA
B	EWALLY-001 EndoWally			EA
C	SCPH-102 Scope Holder			EA
D	30-636 1 PANEL Louvered Panel	35.75"	19.00"	EA
E	30-148 Wall-Mounted Bin Rail	48.00"	3.00"	EA

Plastic Utility Carts

Heavy duty, scratch & stain resistant

Plastic Utility Carts include a lifetime warranty. They will not stain, scratch, dent or rust. They are constructed of strong, heavy-duty, high-density polyethylene. The integral push handle is molded into the top shelf.



Utility Carts

Offered with two or three shelves

Each kit includes two U-handle posts, four casters and either two or three shelves. 500lb. capacity. Casters are 4" diameter (18" shelves) or 5" resilient rubber and each is swivel-style. Shipped as components and quickly assembled with the use of a hammer.



	Item Number	Description	Length	Width	Height	Qty
A	HE-32 BK	2-Shelf Plastic Utility Cart - Black	24.00"	18.00"	33.00"	EA
B	HE-32 GY	2-Shelf Plastic Utility Cart - Grey	24.00"	18.00"	30.00"	EA
C	EU2-1824Z	Utility Cart - Two Shelves	18.00"	24.00"	37.00"	EA
	EU2-1830Z	Utility Cart - Two Shelves	18.00"	30.00"	37.00"	EA
	EU2-1836Z	Utility Cart - Two Shelves	18.00"	36.00"	37.00"	EA
	EU2-2136Z	Utility Cart - Two Shelves	21.00"	36.00"	37.00"	EA
	EU2-2436Z	Utility Cart - Two Shelves	24.00"	36.00"	37.00"	EA
D	EU3-1824Z	Utility Cart - Three Shelves	18.00"	24.00"	37.00"	EA
	EU3-1830Z	Utility Cart - Three Shelves	18.00"	30.00"	37.00"	EA
	EU3-1836Z	Utility Cart - Three Shelves	18.00"	36.00"	37.00"	EA
	EU3-2136Z	Utility Cart - Three Shelves	21.00"	36.00"	37.00"	EA
	EU3-2436Z	Utility Cart - Three Shelves	24.00"	36.00"	37.00"	EA

HushKarts

Quiet mobility & easy to clean design

HushKarts are designed for quiet mobility and easy cleaning with one-piece, double-wall polyethylene construction. Each HushKart offers three shelves with slip resistant surfaces as well as swivel casters with one brake.



	Item Number	Description	Color	Length	Depth	Height	Load Capacity	Qty
A	QCT-225 SB	Open HushKart	Slate Blue	28.00"	16.00"	32.25"	350 lbs.	EA
B	QCT-225 BG	Open HushKart	Beige	28.00"	16.00"	32.25"	350 lbs.	EA
C	QCT-225 BN-FRM	Open HushKart	Brown	28.00"	16.00"	32.25"	350 lbs.	EA
D	QCT-3301 SB	Closed HushKart	Slate Blue	33.125"	20.00"	34.625"	500 lbs.	EA
E	QCT-2301 CG	Open HushKart	Charcoal Gray	33.25"	20.00"	34.625"	500 lbs.	EA

Disposable Cart Covers

Eliminate the expense of laundering covers

Available in clear plastic when you want contents visible and opaque when you want contents out of sight. Made of tough polyethylene plastic, these covers protect carts/contents and eliminate the time and expense of laundering cloth covers.



Item Number	Description	Color	Length	Width	Height	Qty
667-362046-1M RED	Disposable Cart Cover	Red	20.00"	36.00"	46.00"	150
A 668-2 RED	Disposable Cart Cover	Red	24.00"	48.00"	45.00"	50
668-1 TAN	Disposable Cart Cover	Tan	24.00"	36.00"	37.00"	50
B 668-2 TAN	Disposable Cart Cover	Tan	24.00"	48.00"	34.00"	50
668-3 TAN	Disposable Cart Cover	Tan	24.00"	48.00"	48.00"	50
C 667-382648-1M BLUE	Disposable Cart Cover	Blue	38.00"	26.00"	48.00"	50
667-202450-1M	Disposable Cart Cover	Clear	20.00"	24.00"	50.00"	50
666-1 CLEAR	Disposable Cart Cover	Clear	24.00"	27.00"	34.00"	50
668-3 CLEAR	Disposable Cart Cover	Clear	24.00"	48.00"	48.00"	50
667-256161-1M	Disposable Cart Cover	Clear	24.00"	60.00"	60.00"	50
667-252958-1M	Disposable Cart Cover	Clear	25.00"	29.00"	58.00"	50
667-274036-1M	Disposable Cart Cover	Clear	27.00"	40.00"	36.00"	50
667-382648-1M	Disposable Cart Cover	Clear	38.00"	26.00"	48.00"	50

Custom Reusable Cart Covers

Choose size, material, color & closure

Lightweight, flexible and durable custom-made cart covers protect linen from cross contamination during transportation. All covers have reinforced corners and are available with zipper or Velcro tab closures. Both standard and custom-sized rack covers are manufactured with quality materials and workmanship.

Material choices include 200 and 400 denier nylon, bacteria static Staph-Check and clear vinyl. Choose either zipper or Velcro closures and select the sides of the cover for flaps. Select from the many colors available. Custom sizes are also available.



200 & 400 Denier Nylon Color Options:								Staph-Check Color Options:																
Item Number	Details	Width	Length	Height	Qty																			
A SRC-835	Reusable Cart Cover	18.00"	36.00"	54.00"	EA	Red	Royal	Navy	Yellow	White	Kelly	Forest	Black	Vanilla	Powder Blue	Royal Blue	Sand	Black	Pistachio	White	Red	Navy	Forest	Dusty Rose
SRC-836	Reusable Cart Cover	18.00"	36.00"	63.00"	EA																			
SRC-837	Reusable Cart Cover	18.00"	36.00"	74.00"	EA																			
SRC-435	Reusable Cart Cover	24.00"	36.00"	54.00"	EA																			
SRC-436	Reusable Cart Cover	24.00"	36.00"	63.00"	EA																			
SRC-437	Reusable Cart Cover	24.00"	36.00"	74.00"	EA																			
SRC-438	Reusable Cart Cover	24.00"	36.00"	86.00"	EA																			
SRC-440	Reusable Cart Cover	24.00"	48.00"	63.00"	EA																			
SRC-441	Reusable Cart Cover	24.00"	48.00"	72.00"	EA																			
SRC-442	Reusable Cart Cover	24.00"	60.00"	54.00"	EA																			
SRC-443	Reusable Cart Cover	24.00"	60.00"	63.00"	EA																			
SRC-444	Reusable Cart Cover	24.00"	60.00"	72.00"	EA																			



TearGuard Shelf Liner

End ripped wrap due to wire shelving

These liners provide a solid plastic barrier from the sharp edges and bumps of wire shelving and their raised ribs promote good air circulation. They are manufactured of tough polyethylene plastic that can be cleaned with disinfecting solutions and is cart washer safe.



A

Item Number	Details	Length	Width	Qty
A SH1830C	TearGuard Shelf Liner	18.00"	30.00"	EA
SH1836C	TearGuard Shelf Liner	18.00"	36.00"	EA
SH1842C	TearGuard Shelf Liner	18.00"	42.00"	EA
SH1848C	TearGuard Shelf Liner	18.00"	48.00"	EA
SH1860C	TearGuard Shelf Liner	18.00"	60.00"	EA
SH1872C	TearGuard Shelf Liner	18.00"	72.00"	EA
SH2430C	TearGuard Shelf Liner	24.00"	30.00"	EA
SH2436C	TearGuard Shelf Liner	24.00"	36.00"	EA
SH2442C	TearGuard Shelf Liner	24.00"	42.00"	EA
SH2448C	TearGuard Shelf Liner	24.00"	48.00"	EA
SH2460C	TearGuard Shelf Liner	24.00"	60.00"	EA
SH2472C	TearGuard Shelf Liner	24.00"	72.00"	EA



Shelf Mats

Keep items from falling through shelves

Shelf Mats are designed to keep items from falling through open shelving. They are made of durable .035 thick sheets of translucent polyethylene. They can be easily removed and cleaned and available in 12 different sizes, which fit the most commonly used shelf sizes.



A

Item Number	Details	Length	Width	Qty
SM1836CL	Clear Shelf Mat	18.00"	36.00"	5
SM1848CL	Clear Shelf Mat	18.00"	48.00"	5
SM1860CL	Clear Shelf Mat	18.00"	60.00"	5
SM1872CL	Clear Shelf Mat	18.00"	72.00"	5
SM2436CL	Clear Shelf Mat	24.00"	36.00"	5
SM2448CL	Clear Shelf Mat	24.00"	48.00"	5
SM2460CL	Clear Shelf Mat	24.00"	60.00"	5
SM2472CL	Clear Shelf Mat	24.00"	72.00"	5
A SI2136	Shelf Insert	21.00"	36.00"	EA
SI2148	Shelf Insert	21.00"	48.00"	EA
SI2160	Shelf Insert	21.00"	60.00"	EA
SI2172	Shelf Insert	21.00"	72.00"	EA



A Getinge company

Intelligent Solutions
for Instrument Care &
Infection Control

hmark.com

Healthmark, A Getinge company

18600 Malyn Blvd.

Fraser, Michigan 48026

USA

Call Us

Toll Free: 800.521.6224

International: +1.586.774.7600

Fax: 586.491.2113

Email Us

Contact: healthmark@hmark.com

Order Inquiries: hmorders@hmark.com

