




Healthmark Industries Company

LABELING & IDENTITY


A complete line of inventory control, message and customs labels for identification, tracking and communication.



Tray Identification Tags
Page 2




Tray Marking & Tracking
Page 3



Clip-On Basket Tags
Page 4



Label Gun and Labels
Page 8



Color Code Label Tape & Dots
Page 9



Custom Labels
Page 10



Message Labels
Page 12



Custom Signs
Page 15

SURGERY LINE-UP 4-19-08				
O.R.	TIME	SURGEON	STATUS	PROCEDURE
	7:30A	MARKS		
	9:00A	Jones		
	7:30A	STILLS		
	9:30A	HANSON		
	7:00A	KENT		
	8:15A	CARLSON		
	10:00A	KENT		
	8:00A	MASON		
	10:00A	MARKS		

■ O.R. 1 ■ O.R. 3
 ■ O.R. 2 ■ O.R. 4

Insta Grids
Page 20

Call:800-521-6224

Fax:586-491-2113

www.HMARK.com



INSTRUMENT & TRAY TRACKING

Tray Identification Tags

Resin bonded tags designed for tray identification. Tags may be used with bar code labels or hand written notes. All instruments remain with the sterilization basket throughout the collection, decontamination and sterilization process. X-ray detectable.



Code No	Description	Width	Length
RIN001	Stainless Ring, 100 ea.	1.00	
TG001	White, 100 ea.	3.50	2.20
TG011	White, 100 ea.	3.50	1.50
TG012	Tan, 100 ea.	3.50	1.50
TG013	Yellow, 100 ea.	3.50	1.50
TG014	Orange, 100 ea.	3.50	1.50
TG015	Red, 100 ea.	3.50	1.50
TG016	Blue, 100 ea.	3.50	1.50
TG017	Green, 100 ea.	3.50	1.50
TG018	Grey, 100 ea.	3.50	1.50
TG019	Purple, 100 ea.	3.50	1.50
TG020	Black, 100 ea.	3.50	1.50
TG021	Lime, 100 ea.	3.50	1.50
TG022	Pink, 100 ea.	3.50	1.50
TG023	Turquoise, 100 ea.	3.50	1.50
TG024	Mustard, 100 ea.	3.50	1.50
TG025	Terracotta, 100 ea.	3.50	1.50
TG026	Royal Blue, 100 ea.	3.50	1.50
TIE003	Plastic Tie, 100 ea.	0.14	5.91


INSTRUMENT & TRAY TRACKING



Instrument Self Looping Labels



These labels facilitate specific markings and remarks for sterile or contaminated instrument trays. The label is suitable for sterilization and are resistant to washing and disinfection and thus can stay with a tray of instruments all the way through the decontamination process. Packaged 500 per box. Need a special design? Call for a custom label quote.



Code No	Description	Length	Width
403005	White Tray Label	6.50	4.00
403115	White Tray Label	9.50	2.00
403125	Preprinted Tray Label	9.50	2.00
403225	White Tray Label	8.50	2.50
403225 BL	Blue Tray Label	8.50	2.50
403225 GN	Green Tray Label	8.50	2.50
403225 BK	Black Tray Label	8.50	2.50
403225 OR	Orange Tray Label	8.50	2.50
403225 HTK	Label with Month and Dates (500 per box/also available in green, yellow, blue, orange and pink)	8.50	2.50
403225 LNR	"LOANER" White Tray Label	8.50	2.50
403225 PK	Pink Tray Label	8.50	2.50
403225 RD	Red Tray Label	8.50	2.50
403225-RD-BIO	Red Biohazard Label	8.50	2.50
403225 SLTD	White, Slotted for Count Sheet	8.50	2.50
403225 SLTD PK	Pink, Slotted for Count Sheet	8.50	2.50
403225 YL	Yellow Tray Label	8.50	2.50
404225-351-GI	Flexible Endoscope Label	13.00	2.50
STT-001	Self-Tie Tag	11.00	3.75
404225-351-GI	Flexible Endoscope Label	13.00	2.50
EWTC-001	Eye Wash Temperature Check	9.00	5.00

Call:800-521-6224

Fax:586-491-2113

www.HMARK.com




BASKET TAGS



Clip-On Basket Tags

Healthmark's Clip-On Basket Tags (3.5 x 2) are available in 8 different colors to allow users to differentiate instrument sets during reprocessing. Made of tough polypropylene, these tags are able to withstand heat up to 134°C (275 F). Basket tags are compatible with steam, gas, or plasma sterilization. Available in 8 colors: orange, red, blue, yellow, green, black, pink and purple. Basket tags can be purchased in singles or packs of 25.

Code No	Description	Length	Width
CBTOR-1	Clip-On Basket Tag-Orange	3.500	2.000
CBTRD-1	Clip-On Basket Tag-Red	3.500	2.000
CBTBL-1	Clip-On Basket Tag-Blue	3.500	2.000
CBTYL-1	Clip-On Basket Tag-Yellow	3.500	2.000
CBTGN-1	Clip-On Basket Tag-Green	3.500	2.000
CBTBK-1	Clip-On Basket Tag-Black	3.500	2.000
CBTPP-1	Clip-On Basket Tag-Purple	3.500	2.000



Stainless Steel Tray Tags

Available in 3 styles, these stainless steel tags are a great alternative for identifying trays, containers and other key assets. Choose from clip-on, tie-on in round and rectangle.

Code No	Description	Length	Width
01444	Stainless Steel Tag	3.50	1.50
01452	Stainless Steel Badge	0.75	0.00
01454	Clip-on Stainless Steel	3.50	1.50




Laser CAATS

Customize your own basket tags to identify trays, containers and many other items. These custom labels are laser engraved into high quality anodized aluminum. They are available in a variety of colors and sizes.



ROUND TAGS ARE AVAILABLE IN THE FOLLOWING COLORS:



Code No	Color	Length	Width
20RBL	Blue	2.00	2.00
20ROR	Copper Orange	2.00	2.00
20RGN	Green	2.00	2.00
20RPP	Purple	2.00	2.00
20RRD	Red	2.00	2.00
20RSV	Silver	2.00	2.00

Code No	Color	Length	Width
15RBL	Blue	1.50	1.50
15ROR	Copper Orange	1.50	1.50
15RGN	Green	1.50	1.50
15RPP	Purple	1.50	1.50
15RRD	Red	1.50	1.50
15RSV	Silver	1.50	1.50
15RGD	Gold	1.50	1.50



RECTANGULAR TAGS ARE AVAILABLE IN THE FOLLOWING COLORS:



Code No	Color	Length	Width
23RCBL	Blue	2.000	3.000
23RCBLK	Black	2.000	3.000
23RCSV	Silver	2.000	3.000
23RCGD	Gold	2.000	3.000
23RCRD	Red	2.000	3.000



INTERCHANGEABLE & MAGNA TAGS

Interchangeable Tags


Constructed from stainless steel the heat resistant insert holder is manufactured with a durable clip that attaches to various medical trays and mesh wire baskets. The holder is accompanied by a variety of interchangeable color insert tags for easy identification. Made of resilient Polyethersulfone, the interchangeable insert tags are heat resistant up to 275 °F.

Code No.	Color
01498	Stainless Steel Clip-on Tag Holder 2.4" x 1.6"
01948	White Insert Tag
01949	Red Insert Tag
01950	Green Insert Tag
01969	Dark Blue Insert Tag
01951	Light Blue Insert Tag
01952	Orange Insert Tag
01953	Black Insert Tag
01954	Yellow Insert Tag
01990	Pink Insert Tag
01991	Purple Insert Tag



Magna Tags

Quickly label any magnetic surface with the Magna Tags. Attach your own label on the surface of the tag, or use a choose from six different colors and write directly on the erasable surface with a magic marker.



Peel & Stick Insta Grids

Code No.	Color	Width	Height
860	Red	15.33	23.00
863	Green	15.33	23.00
864	Blue	15.33	23.00
865	Black	15.33	23.00

Magnet Tags

Code No.	Color	Width	Height
870	Red	3.50	1.25
871	Orange	3.50	1.25
872	Yellow	3.50	1.25
873	Green	3.50	1.25
874	Blue	3.50	1.25
875	Black	3.50	1.25

SPECIALTY LABELS



Transportation & Quick Turnover Labels



Designed for compliance with OSHA standard CFR 1910.1030 the Transportation Identification Tag is 3.125" x 5.125" with one perforated tab, a green top tab with "CLEAN" in black text, a fluorescent orange/red bottom tab with "DIRTY" in black text, and the removable OSHA approved "Biohazard Label" adhesive backing. The Transportation Identification Tag is produced precisely for transporting materials considered a biohazard, while acting as an essential communication tool in the process.

Code No.	Description
AV-52481	Transportation Label
AV-52482	Transportation Label w/o checklist
HM-52489	Quick Turnover Label

7 Day, 12 Day, 72 hour and 1 Hour Indicator Labels



HTK-7D2

Designed as a visual reminder for healthcare workers, the 7 Day & 12 Day Indicator Hangtime Label features a elapsed time indicator adhered to a 8.5 x 2.5 inch self-looping label. To activate the indicator simply press the blister on the strip, and upon activation a safe dye migrates along a seven or twelve day display run-out window, which shows elapsed time. No preconditioning is required, there are no moving parts, and no power source is required.



HTK-72H



HTK-12D

Operating similar to the 7 & 12 Day Labels, the 72 & 1 Hour indicators will elapse in hour time to represent the amount of time after bedside cleaning that reprocessing must begin.



HTK-1H2



LABEL GUN SYSTEM

Handheld Labeling Systems


Healthmark's Handheld Labeling System has been designed specifically with Central Sterile in mind for the labeling of any items that you need to be autoclaved - from sterile linen packs to instrument trays, and peel pouches. It also has great versatility for labeling in other hospital departments such as dietary and pharmacy.

This system gives you the flexibility to print two lines of information with up to eight characters per line. The available characters include two gregorian dates for printing expiration and autoclave dating with additional bands for load and sterilizer numbers.

Gas-Plasma compatible labels are available. "S" at end of code number is the Gas-Plasma Compatible label.

NEW! Our newest premium label gun prints two-part labels designed to allow the sterilization record label to go with the patient record after use. It also allows you to print more characters on each label.

Event Related




H369023/H369023S

Dual Layer Event Related




15.500.518

Expiration




H369013/H369013S

Dual Layer Expiration



15.500.525

Label Text Colors:




Code No.	Description	Pkg
H-362026	Handheld Label Gun	EACH
H369020	Event related label - red (steam & eto)	CS-12
H369021	Event related label - blue (steam & eto)	CS-12
H369023	Event related label - black (steam & eto)	CS-12
H369024	Event related label - green (steam & eto)	CS-12
H369010	Expiration label - blue (steam & eto)	CS-12
H369012	Expiration label - red (steam & eto)	CS-12
H369013	Expiration label - black (steam & eto)	CS-12
H369014	Expiration label - green (steam & eto)	CS-12
H369020S	Event related sterrad labels - red (sterrad, steam & eto)	CS-12
H369021S	Event related sterrad labels - blue (sterrad, steam & eto)	CS-12
H369023S	Event related sterrad labels - black (sterrad, steam & eto)	CS-12
H369024S	Event related sterrad labels - green (sterrad, steam & eto)	CS-12
H369010S	Expiration labels sterrad - blue (sterrad, steam & eto)	CS-12
H369012S	Expiration labels sterrad - red (sterrad, steam & eto)	CS-12
H369013S	Expiration labels sterrad - black (sterrad, steam & eto)	CS-12
H369014S	Expiration labels sterrad - green (sterrad, steam & eto)	CS-12

Code No.	Description	Pkg
15.400.250	Premium Label Gun	EACH
15.500.515	Event Red Label(steam & eto)	CS-14
15.500.516	Event Green Label (steam & eto)	CS-14
15.500.517	Event Blue Label(steam & eto)	CS-14
15.500.518	Event Black Label(steam & eto)	CS-14
15.500.519	Expiration Red Label(steam & eto)	CS-14
15.500.523	Expiration Green Label(steam & eto)	CS-14
15.500.524	Expiration Blue Label (steam & eto)	CS-14
15.500.525	Expiration Black Label(steam & eto)	CS-14



Color Coded Label Tape

This tape readily accepts ballpoint, pencil or permanent markers. With a service temperature of -70°F to +250°F, it can be used in freezer, refrigerator or incubator. This tape will stay put until you purposely remove it and comes off cleanly with no adhesive transfer. Available in 2 convenient length rolls, 5 widths and 26 colors.



Code No.	Length	Width	Pkg
LT-512	500	0.50	ROLL
LT-534	500	0.75	ROLL
LT-510	500	1.00	ROLL
LT-614	2160	0.25	ROLL
LT-612	2160	0.50	ROLL
LT-610	2160	1.00	ROLL
LT-620	2160	2.00	ROLL
TPEDS-1			1
TPEDS-2			1
TPEDS-3			1

Time Out Sign

Time out! A “Surgical Pause” is the universal protocol for preventing wrong site, wrong procedure, and wrong person surgery. Healthmark’s new stainless steel “Time Out” signs serve as a perfect reminder. The 304 grade stainless is materially compatible with steam, EO & gas-plasma sterilization.



Code No.	Description
SS-TO	Stainless Steel Time Out Sign 5" x 1"

Washer Rack Labels

These labels are used to identify and track washer racks through the washing cycles when being tested with the TOSI™ and other ProFormance™ cleaning verification products. Each label has a unique number that will help determine washer racks that are repeatedly not performing to cleaning standards so the issue can be properly identified. They come in the size 2" x 0.625".



Code No	Description	Pkg
WRL-001	Anodized Aluminum, adhesive label	10



PPE DECALS

PPE Decals




What are they?

Style your Personal Protection Equipment with PPE Decals! These fun designs secure by adhesive directly to your PPE that you are wearing. Whether it is to bring attention to an important cause or showing your pride for your department, these decals are unique and stand out in the crowd. They are very easy to apply, just unpeel and the adhesive back sticks on to the material you are wearing.



Style your PPE with the “Heart of the Hospital” Decals. These stylish decals can be placed directly on the PPE you are already wearing during your shifts.

Pink Ribbon Decals are also available. They are a great way to show support and raise awareness during breast cancer awareness month and Mother’s Day. They soon will be available in other colors.

PPE Decals also come in a variety of styles that can be worn at anytime for fun in the department.



Code No	Length	Width
HRT-33	3.000	3.000
RBN-153	1.500	3.000
BLT-1515BN	15.00	1.53
BWTE-2449PA	2.47	4.94
PCKT-45SG	4.99	4.00
BTTNDWN-3915PK	15.00	3.96
TIE-2416BL	14.66	2.47



CUSTOM VINYL & MAGNETIC LABEL



Custom Vinyl Labels

Create any label specialized for your departments needs. These custom vinyl labels are durable enough to even place into a washer (not steam sterilizable).



Code No.	Description
CVL14	1" x 4"
CVL16	1" x 6"
CVL18	1" x 8"
CVL112	1" x 12"
CVL24	2" x 4"
CVL26	2" x 6"
CVL28	2" x 8"
CVL212	2" x 12"
CVL34	3" x 4"
CVL36	3" x 6"
CVL38	3" x 8"
CVL312	3" x 12"
CVL44	4" x 4"
CVL46	4" x 6"
CVL48	4" x 8"
CVL412	4" x 12"
CVL55	5" x 5"

Code No.	Description
CVL56	5" x 6"
CVL58	5" x 8"
CVL512	5" x 12"
CVL66	6" x 6"
CVL68	6" x 8"
CVL612	6" x 12"
CVL76	7" x 6"
CVL77	7" x 7"
CVL78	7" x 8"
CVL712	7" x 12"
CVL86	8" x 6"
CVL88	8" x 8"
CVL812	8" x 12"
CVL96	9" x 6"
CVL98	9" x 8"
CVL99	9" x 9"
CVL912	9" x 12"

Magnetic Labels

Magnetic labels provide the perfect solution for removing and re-applying labels to metal surfaces. Use them to communicate important messages on cabinets, sinks or any other metal surface in your department.



Code No.	Description
CVL-M14	1" x 4"
CVL-M16	1" x 6"
CVL-M18	1" x 8"
CVL-M112	1" x 12"
CVL-M24	2" x 4"
CVL-M26	2" x 6"
CVL-M28	2" x 8"
CVL-M212	2" x 12"
CVL-M34	3" x 4"
CVL-M36	3" x 6"
CVL-M38	3" x 8"
CVL-M312	3" x 12"
CVL-M44	4" x 4"
CVL-M46	4" x 6"
CVL-M48	4" x 8"
CVL-M412	4" x 12"
CVL-M55	5" x 5"

Code No.	Description
CVL-M56	5" x 6"
CVL-M58	5" x 8"
CVL-M512	5" x 12"
CVL-M66	6" x 6"
CVL-M68	6" x 8"
CVL-M612	6" x 12"
CVL-M76	7" x 6"
CVL-M77	7" x 7"
CVL-M78	7" x 8"
CVL-M712	7" x 12"
CVL-M86	8" x 6"
CVL-M88	8" x 8"
CVL-M812	8" x 12"
CVL-M96	9" x 6"
CVL-M98	9" x 8"
CVL-M99	9" x 9"
CVL-M912	9" x 12"



MESSAGE LABELS

Clean Labels

Effective communication among healthcare professionals is an essential part of success for providing superior care and achieving patient satisfaction, which is why we've created specialty labels. Manufactured to convey key information to healthcare professionals, the clean label is intended to conceal and cover the bio symbol on SST systems when transporting clean medical instruments.



Code No.	Description
CLN-4X4	Clean Label

Multi-purpose Labels

Properly labeling medical equipment can be challenging. The Multi-Purpose Equipment Label was developed as a way to ensure consistency and identify whether medical equipment designated for patient use is either clean, in-use or dirty. The 2" x 3" label has an easy to use pull tab style design with an adhesive backing that doesn't leave residue behind and can be applied to numerous medical devices. The brightly colored label provides high visibility to healthcare professionals, allowing individuals to correctly confirm the status of equipment.



Code No.	Description
HM-CID-613	Multipurpose Equipment Label



High Level Disinfected Label

Designed for containers that hold high-level disinfected or contaminated scopes, the High Level Disinfected Label is 2" x 3" that has a green top tab with "CLEAN" and "HIGH LEVEL DISINFECTED" in black text, as well as a removable self-adhesive orange bottom tab with "BIOHAZARD" in black text that provides space to record the time the scope was pre-cleaned. By labeling containers as "CLEAN" or "BIOHAZARD", the High Level Disinfected Label serves as a communication tool that provides information to healthcare professionals working throughout different departments. The High Level Disinfected Label is available for purchase on a roll of 250.

Code No.	Description
HM-52483-HLD	High Level Disinfected

Preprinted Message Labels

actual size



CSW-xxx



CSW-113



CSW-116



CSW-SSS



CSW-165

PREPRINTED MESSAGE LABELS



Preprinted Message Labels

actual size

**WARNING
GAS STERILIZE
ONLY**

CSW-1212

STERILIZED
Indefinite Shelf Life
Unless Damaged or Opened

CSW-1216

AUTOCLAVED ON _____
EXPIRES ON _____

**RETURN TO CENTRAL SERVICE
FOR RESTERILIZATION**

CSW-1215

**This Tray is
INCOMPLETE**

Missing _____

Date _____

CSW-302

**HIGH LEVEL
DISINFECTED**

Date _____ Initials _____
Process # _____

DC-100214-71-BL

STERILIZED

JANUARY _____ LOAD _____

STERILIZED

FEBRUARY _____ LOAD _____

STERILIZED

MARCH _____ LOAD _____

STERILIZED

APRIL _____ LOAD _____

STERILIZED

MAY _____ LOAD _____

STERILIZED

JUNE _____ LOAD _____

STERILIZED

JULY _____ LOAD _____

STERILIZED

AUGUST _____ LOAD _____

STERILIZED

SEPTEMBER _____ LOAD _____

STERILIZED

OCTOBER _____ LOAD _____

STERILIZED

CSW-xxx
Specify Month

**STERILITY IS GUARANTEED
UNTIL THE PACKAGE IS
DAMAGED OR OPENED**

CSW-254

IMPLANTABLE

HOLD UNTIL _____
DATE _____ TIME _____

CSW-IMP

Code No.	Description	Length	Width	Pkg
CSW-116	EXPIRATION DATE	2.00	0.88	ROLL
CSW-1216	"STERILIZED INDEFINITE SHELF LIFE"	2.50	1.00	RL-1000
CSW-254	STERILITY IS GUARANTEED	2.25	0.88	RL-500
CSW-289	"PLEASE IMMERSER USED INSTRUMENTS IN DETERGENT..."	3.00	0.88	ROLL
CSW-302	THIS TRAY IS INCOMPLETE LABEL - ORANGE	2.25	1.75	ROLL
CSW-IMP	"IMPLANTABLE" LABEL BLACK ON RED	3.00	0.88	RL-320
CSW-MOE	MEDICATION ONLY EXPIRES	2.25	.875	RL-500
CSW-RCS	"RETURN TO CENTRAL SUPPLY" LABEL	2.25	0.88	ROLL
CSW-RRS	REUSABLE RETURN TO SPD	2.25	0.88	ROLL
CSW-SPU	"SINGLE PATIENT USE ONLY- THEN..."	2.25	0.88	ROLL
CSW-SSS	STERRAD SYSTEM STERILIZED	2.25	0.88	ROLL
DC-100214-71-BL	HIGH LEVEL DISINFECTED	2.25	1.00	ROLL

Call:800-521-6224

Fax:586-491-2113

www.HMARK.com



CUSTOM LABELS



CUSTOM LABEL FORM

DESIGN YOUR OWN LABEL BY COMPLETING THE STEPS BELOW

1: Adhesive Type: Permanent Removable Ultra Removable

Please note that adhesives may not stick to porous materials, such as wrap.

Custom labels require a minimum order of 2,000





2: Sterilization Compatibility: Steam EO, EtO Gas Plasma

3: Select a material and text color from the options below. Material Color: _____ Text Color: _____

1: WHITE	2: YELLOW	3: GREEN	4: RED	5: ORANGE	6: BLUE	7: PINK	8: TAN	9: SALMON
10: COPPER	11: ROSE	12: VIOLET	13: LAVENDER	14: GRAY	15: LIME	16: AQUA	17: CHARTREUSE	18: GOLD
19: SILVER	20: MAGENTA	21: DARK RED	22: MAROON	23: DARK BROWN	24: DARK BLUE	25: DARK GREEN	26: BLACK	

4: Select a label size by sketching your layout on the desired size below.

Interested in a size not shown? Please specify the size: _____ x _____ and explain your label below:



5: Please provide your information:

Name: _____

Job Title: _____

Department: _____

Facility: _____

Facility Address: _____

Phone Number: _____

Email: _____

Signature: _____ Date: _____

FAX YOUR COMPLETED FORM TO **586-491-2113** AND A CUSTOMER SERVICE REPRESENTATIVE WILL CONTACT YOU ABOUT YOUR ORDER.

FLOOR SIGNS



Floor Signs

Floor signs are designed to communicate important information to employees and visitors of your facility. These new no-lip profile, UV protected industrial laminate materials, these hi-impact signs will last up to 7 years without losing their top-layer built in traction or their color.



FS-SLW



FS-PPEREQB



FS-SLIP



FS-FA



FS-FE-DNB



FS-EW



FS-KAC



FS-NOTCDNE



FS-RECPAPR



FS-EXWY



FS-FLL



Floor Signs



FS-RECPLAS



FS-RETRNSTK



FS-TRSH-CN



FS-STP-AUTH



FS-CAUTSTP



FS-STP-BAS



STC01



STF01



STT01



STA01



ST4STRP01

FINE SIGNS



FINE SIGNS

Select from our stock signs or design your own (see next page) to alert staff and/or visitors to important information.



FS1624-TO




FS2412-WC



FS1624-MSK



FS2424-FLAM



FS1224-AE



FS2406-CATH



FS1624-BIO



FS2412-EW



FS2406-P&P



FS1224-LIF



FS2424-KAC



FS2424-STP



FS2412-CEIL



FS2412-PPE



FS2424-RSTR



FS2424-LABCOAT



FINE SIGNS

Design Your Own Sign

With a few easy steps, you can design a sign of your own.

- Step 1 - Choose from these standard sign sizes (call for custom sizes).
- Step 2 - Select from 6 Background (poster) Colors.
- Step 3 - Select from 6 Foreground (vinyl) Colors.
- Step 5 - Sketch your layout within the appropriate guides for you size poster (below).
- Step 6 - Fax this form to 586-774-6473. Indicate the quantity desired and voila, you have your own sign!

1 SIZE

Select Width:


- 6"
- 12"
- 16"
- 24"

Select Height:


- 6"
- 12"
- 16"
- 24"

FS0606	6" wide x 6" high	FS1606	16" wide x 6" high
FS0612	6" wide x 12" high	FS1612	16" wide x 12" high
FS0616	6" wide x 16" high	FS1616	16" wide x 16" high
FS0624	6" wide x 24" high	FS1624	16" wide x 24" high
FS1206	12" wide x 6" high	FS2406	24" wide x 6" high
FS1212	12" wide x 12" high	FS2412	24" wide x 12" high
FS1216	12" wide x 16" high	FS2416	24" wide x 16" high
FS1224	12" wide x 24" high	FS2424	24" wide x 24" high


3 BACKGROUND COLOR



4 FOREGROUND COLOR



5 LAYOUT



Safe-T-Cals

SAFE-T-CALS


These easy to apply/remove Safe-T-Cals are a fun and effective way to communicate messages throughout your department. They can be made in a variety of sizes. Choose your own art and message!



DESIGN YOUR OWN WALL GRAPHIC SIGN

With a few easy steps, you can design a sign of your own.




- Step 1 - Choose a high resolution image and size
- Step 2 - Select your own quote bubble comments
- Step 3 - Talk to a Customer Service Representative for a custom quote for your wall graphic



Code No.	Description
WG-PPESEC	PPE Security Sign
WG-WAI	Smiley Face Decontam Reminder
WG-PPEPHARM	PPE Hazardous Drugs



INSTA GRIDS



INSTA GRIDS

Insta Grids are grids that are ready to use in seconds! They are available for use on a magnetic board or peel & stick version for use on any whiteboard. Just stick on the board creating an instant grid. Each square is perfectly shaped and presents a professional appearance. Insta Grids can be customized to any application including creating larger rectangles by cutting grid cells to size. Connect grids together to add columns and rows as needed. Insta Grids are available in six colors. Consider ordering different colors to color code your data. One kit contains six 15" x 23" grid panels to complete a 4' x 4' grid board.

SURGERY LINE-UP				
O.R.	TIME	SURGEON	STATUS	4-19-08 PROCEDURE
█	7:30A	MARKS		
█	9:00A	Jones		
█	7:30A	STILLS		
█	9:30A	HANSON		
█	7:00A	KENT		
█	8:15A	CARLSON		
█	10:00A	KENT		
█	8:00A	MASON		
█	10:00A	MARKS		

█ O.R. 1 █ O.R. 3
 █ O.R. 2 █ O.R. 4

MAGNETIC INSTA GRIDS

Magnetic board or peel & stick version for use on any whiteboard.

Code No.	Color	Width	Height
850	Red	15.33	23.00
851	Orange	15.33	23.00
852	Yellow	15.33	23.00
853	Green	15.33	23.00
854	Blue	15.33	23.00
855	Black	15.33	23.00