



## **FOR IMMEDIATE RELEASE:**

### **CONTACTS:**

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Healthmark Offers New Humipak for Robotic Instruments**

Fraser, MI August 29, 2024 – Healthmark, A Getinge Company has announced they are now offering the Humipak for Robotic Instruments to its SST Instrument Retrieval product line.

Manufactured from polyester-polypropylene films with a highly hydrophilic, blue viscose inner liner, the Humipak for Robotic Instruments consists of a layer of highly absorbent material sandwiched between two layers of waterproof film and are offered in the following sizes: 8.3 x 32 In. and 16.5 x 32 In. Designed as single use, the 8.3 x 32 In. Humipak is designed for a single robotic instrument, whereas the 16.5 x 32 In. Humipak is for four robotic instruments. Simply, place the instruments inside the Humipak, add the specified amount of water to the absorbent layer, and seal with the peel away adhesive strip. This creates a watertight, moist atmosphere that will prevent organics from drying over an extended period. The Humipak for Robotic Instruments are available for purchase in a pack of twenty-five and a case of one hundred.

### **About Healthmark, A Getinge Company**

*Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries' mission is to continue to innovate, continue to support, and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.*