



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Healthmark, A Getinge company offers New Angled Toothbrush Style Brush

October 31, 2024 Fraser, MI - Healthmark, A Getinge company has announced that they have added the Angled Toothbrush Style Brush to its ProSys™ Instrument Care product line.

This is a general cleaning brush that can be used to clean the external surface of the distal end and forceps elevator mechanism of endoscopes as well as the external surfaces and openings of camera heads. The overall style of the toothbrush allows it to be used for cleaning multiple types of instruments and devices in decontamination departments, during manual cleaning.

Designed as single-use, the handle is angled at the bristle head adding extra manual pressure. The two different bristle heights assist in cleaning uneven surfaces and allow the taller, soft bristles at the top edge to clean delicate areas. Made up of a polypropylene handle and nylon bristles, the Angled Toothbrush Style Brush is available for purchase in a pack of ten and a case of one hundred.

Go to hmark.com for more information or call 1-800-521-6224 today!

About Healthmark, A Getinge company

Since 1969, Healthmark has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries' mission is to continue to innovate, continue to support, and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.